# **PCT**

#### WORLD INTELLECTUAL PROPERTY ORGANIZATION International Bureau



#### INTERNATIONAL APPLICATION PUBLISHED UNDER THE PATENT COOPERATION TREATY (PCT)

(51) International Patent Classification <sup>5</sup> :		(11) International Publication Number: WO 94/20539
C07K 15/00, C12P 21/08, G01N 33/577, 33/68	A1	(43) International Publication Date: 15 September 1994 (15.09.94)
(21) International Application Number: PCT/US (22) International Filing Date: 4 March 1994 ( (30) Priority Data: 027,070 4 March 1993 (04.03.93) (71) Applicant: CREATIVE BIOMOLECULES, INC. [US South Street, Hopkinton, MA 01748 (US). (72) Inventors: RUEGER, David, C.; 19 Downey Street, ton, MA 01748 (US). JONES, William, K.; 35 S Street, Brookline, MA 02146 (US). TUCKER, R 132 Robert Road, Holliston, MA 01746 (US). MANN, Hermann; 25 Summer Hill Road, Medv 02053 (US). OZKAYNAK, Engin; 44 Purdue Dr ford, MA 01757 (US). SAMPATH, Kuber, T.; Street, Meday, MA 02053 (US). (74) Agent: KELLEY, Robin, D.; Patent Dept., BioMolecules, Inc., 45 South Street, Hopkinton, M (US).	04.03.9 US/US]; 4 Hopkin Saint Pa obert, F OPPEI way, Mi 6 Sprin Creativ	HU, JP, KP, KR, KZ, LK, LV, MG, MN, MW, NO, NZ, PL, RO, RU, SD, SK, UA, UZ, VN, European patent (AT, BE, CH, DE, DK, ES, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE), OAPI patent (BF, BJ, CF, CG, CI, CM, GA, GN, MIL, MR, NE, SN, TD, TG).  Published  With international search report.  Before the expiration of the time limit for amending the claims and to be republished in the event of the receipt of amendments.  R. A.

#### (54) Title: METHOD AND COMPOSITIONS FOR RECOMBINANT OSTEOGENIC PROTEIN PRODUCTION

### (57) Abstract

Disclosed are novel compositions of osteogenic proteins constituting soluble forms of these proteins, and methods and compositions for distinguishing between soluble and mature forms of these proteins.

# FOR THE PURPOSES OF INFORMATION ONLY

Codes used to identify States party to the PCT on the front pages of pamphlets publishing international applications under the PCT.

	4		** · · · · · * **		A Processing Control of the Control
AT	Austria	GB	United Kingdom	MIR	Mauritania
ΑÜ	Australia	GE	Georgia	MW	Malawi
BB	Barbados	GN	Guinea	NE	Niger
BE	Belgium	GR	Greece	NL	Netherlands
BF	Burkina Faso	HU	Hungary	NO	Norway
BG	Bulgaria	Œ	Ireland	NZ	New Zealand
BJ	Benin	IT	Italy .	PL	Poland
BR	Brazil	JP	Japan	PT	Portugal
BY	Belarus	KE	Kcnya	RO	Romania
CA	Canada	KG	Kyrgystan	RU	Russian Federation
CF	Central African Republic	KP	Democratic People's Republic	SD	Sudan
CG	Congo		of Korea	SE	Sweden
CH	Switzerland	KR	Republic of Korea	SI	Slovenia
CI	Côte d'Ivoire	KZ	Kazakhstan	SK	Slovakia
CM	Cameroon	LI	Liechtenstein	SN	Senegal
CN	China	LK	Sri Lanka	TD	Chad
CS	Czechoslovakia	LU	Luxembourg	TG	Togo
CZ	Czech Republic	LV	Latvia	TJ	Tajikistan
DE	Germany	MC	Monaco	TT	Trinidad and Tobago
DK	Denmark	MD	Republic of Moldova	UA	Ukraine
ES	Spain	MG	Madagascar	US	United States of America
FI	Pinland	ML	Mali	UZ	Uzbekistan
FR	France	MIN	Mongolia	VN	Vict Nam
GA	Gabon				

-1-

#### "METHODS AND COMPOSITIONS FOR RECOMBINANT OSTEOGENIC PROTEIN PRODUCTION"

#### Relation to Related Applications

This application is a continuation-in-part of
USSN 08/027,070, filed March 4, 1993, which is a
continuation-in-part of USSN 07/841,646, filed February
21, 1992, now US Patent No. 5,266,683, the disclosures
of which are incorporated herein by reference.

#### 10 Field of the Invention

The present invention relates generally to osteogenic proteins and, more particularly, to methods and compositions for their production and purification.

15

#### Background of the Invention

Osteogenic proteins are well known and described in the art. See, for example, U.S. Pat. Nos. 4, 968,590; 5,011,691; 5,018,753 and 5,266,683, as well as various scholarly articles published in the scientific literature. See, for example, Ozkaynak et al. (1990) EMBO J 9:2085-2093; Ozkaynak et al., J.Biol. Chem. 267:13198-13205; Sampath et al. (1993) PNAS 90: 6004-6008; Wozney et al. (1988) Science 242: 1528-1534; Wang et al. (1988) PNAS 85:9484-9488; Wang et al. (1990) PNAS 87:2220-2224, and Celeste, et al. (1990) PNAS 87:9843-9847. The art has described how to isolate osteogenic proteins from bone and how to identify genes

- 2 -

encoding these protein and express them using recombinant DNA technology.

The proteins that define the class of true 5 osteogenic proteins constitute a group of proteins sharing a number of conserved structural characteristics. Each protein, on its own, can induce endochondral bone formation in a mammal when properly folded, dimerized and disulfide bonded to produce a dimeric species having the appropriate three 10 dimensional conformation, and without requiring the addition of other osteogenic or non-osteogenic proteins. Typically, osteogenic proteins are provided to a site for bone induction in a mammal in association 15 with a suitable matrix having the appropriate conformation to allow the infiltration, proliferation and differentiation of migrating progenitor cells. The construct of osteogenic protein adsorbed to the surfaces of a suitable matrix is generally referred to 20 in the art as an osteogenic device. The proteins can be isolated from bone or, preferably, the gene encoding the protein is produced recombinantly in a suitable host cell. Methods for the production of osteogenic proteins and formulations of osteogenic devices are described in detail in the art. See, for example, 25 U.S.Pat. Nos. 5,011,691, or 5,266,683, the disclosures of which are incorporated hereinabove by reference.

Improved methods for the recombinant expression of osteogenic proteins is an ongoing effort in the art. It is an object of this invention to provide an improvement in the methods for producing and purifying osteogenic proteins having high specific activity, and for formulating osteogenic devices comprising these

- 3 -

proteins. Still another object is to provide means for distinguishing between the soluble form of the protein and the mature osteogenic protein species typically utilized in formulating osteogenic devices, and to 5 provide polyclonal and monoclonal antibodies capable of distinguishing between these various species. Another object is to provide methods for producing antibodies which can recognize both forms of the protein. another object is to provide methods for monitoring each and all these forms of the protein in a fluid, including serum and production media. U.S. Patent No. 4,857,956 and Urist et al. (1984) Proc. Exp. Bio. Med. 176: 472-475, describe a serum assay for detecting a protein purported to have osteogenic activity. protein is not a member of the family of osteogenic proteins described herein, differing in molecular weight, structural characteristics and solubility from these proteins.

#### 20 Summary of the Invention

It has now been discovered that osteogenic protein as defined herein below, when secreted into cultured medium from mammalian cells contains as a significant fraction of the secreted protein a soluble form of the protein, and that this soluble form comprises the mature dimeric species, including truncated forms thereof, noncovalently associated with at least one, and preferably two pro domains. It further has been discovered that binding partners having specific binding affinity for an epitope on an osteogenic protein or precursor polypeptide chain can be used to discriminate between these two forms of the protein. Preferably the binding partner is a protein. In one

- 4 -

preferred embodiment the binding protein is an antibody, which may be monoclonal, polyclonal or biosynthetically produced. These binding partners may be used as part of a purification scheme to selectively 5 isolate the mature or the soluble form of the protein. as well as to quantitate the amount of mature and soluble forms produced. The antibodies may be used as part of a production protocol to monitor the pharmacological purity of an osteogenic protein 10 preparation for therapeutic or other clinical applications. Specifically, a method now is provided herein to ensure that only one desired form of a protein is present in a composition. In addition, binding partners can be produced which recognize both 15 protein forms and which may be used to advantage to monitor the quantity of total protein in solution. These binding partners also may be used as part of diagnostic treatments to monitor the concentration of osteogenic protein in solution in a body and to detect 20 fluctuations in the concentration of the proteins in their various forms.

The foregoing and other objects, features and advantages of the present invention will be made more apparent from the following detailed description of the invention.

## Brief Description of the Drawings

30 Fig. 1 is a schematic representation of an osteogenic protein polypeptide chain as expressed from a nucleic acid encoding the sequence, wherein the cross-hatched region represents the signal sequence; the stippled region represents the pro domain; the

- 5 -

hatched region represents the N-terminus ("N-terminal extension") of the mature protein sequence; and the open region represents the C-terminal region of the mature protein sequence defining the conserved seven cysteine domain, the conserved cysteines being indicated by vertical hatched lines;

Fig.2 lists the sequences of the N-terminal extensions of the mature forms of various osteogenic 10 proteins; and

Fig. 3 is a gel filtration column elution profile of a soluble osteogenic protein (OP-1) produced and purified from a mammalian cell culture by IMAC,

15 S-Sepharose and S-200HR chromatography in TBS (Tris-buffered saline), wherein V<sub>O</sub> is the void volume, ADH is alcohol dehydrogenase (MW 150 kDa), BSA is bovine serum albumin (MW 67 kDa), CA is carbonic anhydrase (MW 29kDa) and CytC is cytochrome C (MW 12.5 kDa).

#### Detailed Description

A soluble form of true osteogenic proteins now has
been discovered wherein the protein consists
essentially of the amino acid sequence of the protein.
The soluble form is a non-covalently associated complex
comprising the pro domain or a fragment thereof,
noncovalently associated or complexed with a dimeric
protein species having osteogenic activity, each
polypeptide of the dimer having less than 200 amino
acids and comprising at least the C-terminal six, and
preferably seven cysteine skeleton defined by residues
335-431 and 330-431, respectively, of Seq. ID No. 1.

- 6 -

Preferably, the polypeptide chains of the dimeric species comprise the mature forms of these sequences, or truncated forms thereof. Preferred truncated forms comprise the intact C-terminal domain and at least 10 amino acids of the N- terminal extension sequence e.g., preferably at least the sequence defined by residues 320-330 of Seq. ID 1. The soluble forms of these osteogenic proteins may be isolated from cultured cell medium, a mammalian body fluid, or may be formulated in vitro.

In vivo, under physiological conditions, the prodomain may serve to enhance the transportability of the proteins, and/or to protect the proteins from proteases and scavenger molecules, including antibodies. The prodomains also may aid in targeting the proteins to tissue, e.g., to bone.

The proteins contemplated by the invention and referred to herein as "osteogenic proteins" are true osteogenic proteins capable, on their own, of inducing endochondral bone formation when implanted in a mammal in association with a matrix, without requiring the addition of other osteogenic or non-osteogenic proteins. A detailed description of these proteins appears, for example, in U.S. Pat. Nos. 4968,590, 5,011,691 and USSN 841,646, and includes references to various members of the protein family identified to date. These family members include OP1, OP2, and the proteins referred to in the art as "bone morphogenic proteins": BMP-2, BMP-3, BMP-4, BMP-5, BMP-6 and BMP-9, as well as various known species variants, including Vgr, Vgl, 60A and DPP, and biosynthetic

- 7 <del>-</del>

osteogenic constructs, including COP1, 3, 5, 7 and 16.

The members of this family of proteins, which are a subclass of the TGF-\$\beta\$ super-family of proteins, share 5 characteristic structural features, represented schematically in Fig. 1, as well as substantial amino acid sequence homology in their C-terminal domains, including a conserved seven cysteine structure. As illustrated in the figure, the proteins are translated as a precursor polypeptide sequence 10, having an N-10 terminal signal peptide sequence 12, (the "pre pro" region, indicated in the figure by cross-hatching), typically less than about 30 residues, followed by a "pro" region 14, indicated in the figure by stippling, 15 and which is cleaved to yield the mature sequence 16. The mature sequence comprises both the conserved Cterminal seven cysteine domain 20, and an N-terminal sequence 18, referred to herein as an N-terminal extension, and which varies significantly in sequence 20 between the various osteogenic proteins. Cysteines are represented in the figure by vertical hatched lines 22. The polypeptide chains dimerize and these dimers typically are stabilized by at least one interchain disulfide bond linking the two polypeptide chain subunits. The mature subunits produced from mammalian 25 cells typically have molecular weights in the range of about 15-23 kD, depending on the degree of glycosylation and N-terminal truncation. The dimeric species then typically have a molecular weight in the range of about 30-40 kD. 30

The signal peptide is cleaved during or soon after translation, at a cleavage site that can be predicted

- 8 -

in a given sequence using the method of Von Heijne
((1986) Nucleic Acids Research 14:4683-4691.) The
"pro" form of the protein subunit, 24, in Fig. 1,
includes both the pro domain and the mature domain,

peptide bonded together. Typically, this pro form is
cleaved while the protein is still within the cell, and
the pro domain remains noncovalently associated with
the mature form of the subunit to form a soluble
species that appears to be the primary form secreted
from cultured mammalian cells. Typically, previous
purification techniques utilized denaturing conditions
that disassociated the complex.

Other soluble forms of osteogenic proteins secreted

from mammalian cells include dimers of the pro forms of these proteins, wherein the pro domain is not cleaved from the mature domain, and "hemi-dimers", wherein one subunit comprises a pro form of the polypeptide chain subunit and the other subunit comprises the cleaved

mature form of the polypeptide chain subunit (including truncated forms thereof), preferably noncovalently associated with a cleaved pro domain.

The isolated pro domain typically has regions of

25 hydrophobicity, as determined both by analysis of the
sequence and by characterization of its properties in
solution. The isolated pro domains alone typically are
not fully soluble in aqueous solutions. Accordingly,
without being limited to any given theory, the non
30 covalent association of the cleaved pro domains with
the mature osteogenic protein dimeric species may
involve interaction of a hydrophobic portion of a given
pro domain with a corresponding hydrophobic region on
the dimeric species, the interaction of which

- 9 -

effectively protects or "hides" an otherwise exposed hydrophobic region of the mature dimer from exposure to aqueous environments, enhancing the affinity of the mature dimer species for aqueous solutions.

5

Like the osteogenic proteins described herein,  $TGF-\beta$  also has a pro domain which associates non-covalently with the mature  $TGF-\beta$  protein form. However, unlike the osteogenic proteins described herein, the  $TGF-\beta$  pro domain contains numerous cysteines and forms disulfide bonds with a specific binding protein. The  $TGF-\beta 1$  pro domain also is phosphorylated at one or more mannose residues, while the osteogenic protein pro domains typically are not.

15

As described above, the active form of osteogenic protein, exemplified herein below by OP-1, is known to comprise a dimeric species composed of the mature sequence (e.g., amino acids 293-431 of Seq. ID NO. 1) 20 or a truncated form thereof, appropriately disulfide bonded to produce an osteogenic dimeric species. These osteogenic proteins, in their mature forms, are neutral to basic proteins (e.g., pI in the range of about 7 to 8) and are, to varying degrees, relatively insoluble 25 under physiological conditions. It now has been discovered that these proteins are secreted from mammalian cells in a soluble form and that this form comprises the mature dimeric species (also referred to herein as the "purified species") associated 30 noncovalently with one or more copies of the pro This form of the protein likely is the form which is present in the bloodstream. Unlike the latent form of TGF-β, this form of osteogenic protein does not inhibit activity. The soluble form itself may be

- 10 -

active or, it may dissociate to release the mature dimeric species when the protein has arrived at a target tissue, such as bone tissue.

5 Thus antibodies now may be generated which recognize an osteogenic protein of interest, and these antibodies then used to monitor culture medium or endogenous levels of osteogenic protein in a body fluid, such as serum, whole blood or peritoneal fluid. 10 Preferably, the antibody has binding specificity for the soluble form. Such antibodies may be generated by using the pro domain or a portion thereof as the antigen, or, preferably, the soluble complex itself. The pro domain may be preferably obtained by isolating 15 the soluble complex and then separating the noncovalently associated pro domain from the mature domain using standard procedures, e.g., by separating the complex components by chromatographic means, preferably by ion-exchange chromatography in the 20 presence of a denaturant, e.g., 6M urea. Alternatively, the pro protein in its monomeric form may be used as the antigen and the candidate antibodies screened by western blot or other standard immunoassay for those which recognize the pro form or soluble form 25 of the protein of interest, but not the mature form. Where antibody capable of identifying both soluble and mature forms of the protein is desired, the complex itself preferably is used as the antigen source. Details of antibody production and exemplary 30 immunoassays are provided below. Also provided is an example for detecting soluble osteogenic protein in a body fluid sample.

- 11 -

The proteins contemplated to be useful in the methods and compositions of the invention include forms having varying glycosylation patterns and varying N-termini. They may be naturally occurring or 5 biosynthetically derived, and may be produced by expression of recombinant DNA in procaryotic or eucaryotic host cells. The proteins are active as a single species (e.g., as homodimers), or combined as a mixed species. Useful sequences and eucaryotic and procaryotic expression systems are well described in the art. See, for example, US Patent Nos. 5,061,911 and 5,266,683 for useful expression systems. Useful sequences are recited in US Pat Nos. 4, 968,590; 5,011,691; 5,018,753 and 5,266,683, and Ozkaynak et al. 15 (1990) EMBO J 9:2085-2093; Ozkaynak et al., J.Biol. Chem. 267:13198-13205; Sampath et al. (1993) PNAS 90: 6004-6008; Wozney et al. (1988) Science 242: 1528-1534; Wang et al. (1988) PNAS 85:9484-9488; Wang et al. (1990) PNAS 87:2220-2224, Celeste, et al. (1990) PNAS 20 87:9843-9847 Weeks et al. (1987) Cell 51:861-867; Padgett, et al. (1987) Nature 325 81-84; Wharton et al. (1991) PNAS 88:9214-9218; Lyons et al. (1989) PNAS 86:4554-4558 and PCT international application WO93/00432 for OP1, OP2, DPP, 60A, Vg1, Vgr-1 and the 25 BMP-2-6, and BMP -9 proteins. Accordingly, these proteins, including allelic, species and other naturally-occurring and biosynthetic sequence variants thereof are contemplated to be useful in the inteant application. Other useful sequences include 30 biosynthetic constructs including, without limitation any of the sequences referred to in U.S. Patent No. 5,011,691 as COP-1, -3, -5, -7, -16; and chimeric constructs created by combining sequences from two or more different osteogenic proteins. As will be

domain has greater than 60% identity, and preferably

domain has greater than 60% identity, and preferably greater than 65% identity with the amino acid sequence of OPS (e.g., residues 335-431 of Seq. ID No. 1).

In another preferred aspect, the invention contemplates osteogenic proteins comprising species of polypeptide chains having the generic amino acid sequence herein referred to as "OPX" which accommodates the homologies between the various identified species of the osteogenic OP1 and OP2 proteins, and which is described by the amino acid sequence presented below and in Sequence ID No. 5.

	Cys	Xaa	Xaa	His	Glu	Leu	Tyr	Val	Ser	Phe
25	1				5					10
	Xaa	Asp	Leu	Gly	Trp	Xaa	Asp	Trp	Xaa	Ile
					15					20
	Ala	Pro	Xaa	Gly	Tyr	Xaa	Ala	Tyr	Tyr	Cys
					25					30
30	Glu	Gly	Glu	Cys	Xaa	Phe	Pro	Leu	Xaa	Ser
					35					40
	Xaa	Met	Asn	Ala	Thr	Asn	His	Ala	Ile	Xaa
					45					50

- 13 -

Gln Xaa Leu Val His Xaa Xaa Xaa Pro Xaa 55 60 Xaa Val Pro Lys Xaa Cys Cys Ala Pro Thr 65 70 5 Xaa Leu Xaa Ala Xaa Ser Val Leu Tvr Xaa 75 80 Asp Xaa Ser Xaa Asn Val Ile Leu Xaa Lys 85 90 Xaa Arg Asn Met Val Val Xaa Ala Cys Gly 10 95 100 Cys His,

and wherein Xaa at res. 2 = (Lys or Arg); Xaa at res. 3 = (Lys or Arg); Xaa at res. 11 = (Arg or Gln); 15 Xaa at res. 16 = (Gln or Leu); Xaa at res. 19 = (Ile or Val); Xaa at res. 23 = (Glu or Gln); Xaa at res. 26 = (Ala or Ser); Xaa at res. 35 = (Ala or Ser); Xaa at res. 39 = (Asn or Asp); Xaa at res. 41 = (Tyr or Cys); Xaa at res. 50 = (Val or Leu); Xaa at res. 52 = (Ser or 20 Thr); Xaa at res. 56 = (Phe or Leu); Xaa at res. 57 = (Ile or Met); Xaa at res. 58 = (Asn or Lys); Xaa at res. 60 = (Glu, Asp or Asn); Xaa at res. 61 = (Thr, Ala or Val); Xaa at res. 65 = (Pro or Ala); Xaa at res. 71 = (Gln or Lys); Xaa at res. 73 = (Asn or Ser); Xaa at 25 res. 75 = (Ile or Thr); Xaa at res. 80 = (Phe or Tyr); Xaa at res. 82 = (Asp or Ser); Xaa at res. 84 = (Ser or Asn); Xaa at res. 89 = (Lys or Arg); Xaa at res. 91 = (Tyr or His); and Xaa at res. 97 = (Arg or Lys).

In still another preferred aspect, the invention contemplates osteogenic proteins encoded by nucleic acids which hybridize to DNA or RNA sequences encoding the C-terminal seven cysteine domain of OP1 or OP2 under stringent hybridization conditions. As used

- 14 -

herein, stringent hybridization conditions are defined as hybridization in 40% formamide, 5 X SSPE, 5 X Denhardt's Solution, and 0.1% SDS at 37°C overnight, and washing in 0.1 X SSPE, 0.1% SDS at 50°C. (See, for example, Molecular Cloning: A Laboratory Manual, Maniatis et al., eds. 2d.ed., Cold Spring Harbor Press, Cold Spring Harbor, 1989.)

Useful pro domains include the full length pro domains described below, as well as various truncated 10 forms hereof, particularly truncated forms cleaved at proteolytic Arg-Xaa-Xaa-Arg cleavage sites. example, in OP-1, possible pro sequences include sequences defined by residues 30-292 (full length form); 48-292; and 158-292. Soluble OP-1 complex 15 stability is enhanced when the pro domain comprises the full length form rather than a truncated form, such as the 48-292 truncated form, in that residues 30-47 show sequence homology to the N-terminal portions of other 20 osteogenic proteins, and are believed to have particular utility in enhancing complex stability for all osteongeic proteins. Accordingly, currently preferred pro sequences are those encoding the full length form of the pro domain for a given morphogen 25 (see below). Other pro sequences contemplated to have utility include biosynthetic pro sequences, particularly those that incorporate a sequence derived from the N-terminal portion of one or more osteogenic protein pro sequences.

30

A brief description of OP-1 is described below, followed by examples disclosing how to isolate soluble forms of these proteins and how to generate and identify antibodies having specificity for the mature

- 15 -

or "purified" form, and the soluble or "media" form, or for both forms of proteins. In a particularly preferred embodiment, antibodies, or other binding proteins which recognize the soluble form of the 5 protein do not also recognize the precursor form of the pro domain peptide alone. As will be appreciated by those having ordinary skill in the art, with this disclosure it is possible to more accurately identify and/or quantitate the amount of a given, preferred form 10 of a recombinantly produced osteogenic protein present in the culture media, than by previously available methods, e.g., methods which relied on antibody specific for epitopes present only on the mature form or soluble form of the protein. It also now is 15 possible to accurately isolate a desired form of the protein. For example, one can preferentially isolate the soluble complex form by passing the culture media over an affinity column that has bound antibody with binding specificity only for the pro domain form, and 20 then selectively desorbing the bound protein using standard procedures to modify binding conditions. thereby allowing selective isolation of the complex. For example, one can modify binding conditions by using low pH, denaturants, or by competition with a peptide specific for the antibody binding site. It also is 25 anticipated that the antibodies and protocols can be used to identify both forms of the protein in a solution. A particularly useful application of the invention is as part of a protocol for monitoring the 30 pharmacological purity of an osteogenic protein composition to be used for clinical applications.

While the examples demonstrate the utility of the invention by means of an illustrative protein, OP1, it

is intended to be understood that the methods and compositions taught herein can be extended without undue experimentation to other members of the family of osteogenic proteins. Similarly, while the examples are directed to antibodies as the binding partner having specificity for an epitope on an osteogenic protein, and immunoassays as the detection protocols, any binding partner, particularly any binding protein capable of providing the same discriminatory power as 10 the antibodies described herein are contemplated. Moreover, while only monoclonal and polyclonal antibodies are described in detail, other antibody forms, including biosynthetics such as single chain constructs, referred to in the art as "sFv's" also are 15 contemplated to be within the scope of the invention.

OP1 - Refers generically to the family of osteogenically active proteins produced by expression of part or all of the hOP1 gene. Also referred to in related applications as "OPI" and "OP-1".

hOP1-PP - Amino acid sequence of human OP1 protein (prepro form), Seq. ID No. 1, residues 1-431. Also referred to in related applications as "OP1-PP" and "OPP".

OP1-18Ser - Amino acid sequence of mature human OP1 protein, Seq. ID No. 1, residues 293-431. N-terminal amino acid is serine. Originally identified as

30 migrating at 18 kDa on SDS-PAGE in COS cells.

Depending on protein glycosylation pattern in different host cells, also migrates at 23kDa, 19kDa and 17kDa on

- 17 -

SDS-PAGE. Also referred to in related applications as "OP1-18."

OPS - Amino acid sequence defining the C-terminal six 5 cysteine domain, residues 335-431 of Seq. ID No. 1.

OP7 - Amino acid equence defining the C-terminal seven cysteine domain, residues 330-431 of Seq. ID No. 1.

- Physical and Antigenic Structure of Soluble OP-1 and Mature, Dimeric OP-1
  - la. Solubility of Mature, Dimeric OP-1
- 15 The solubility properties of purified mature OP-1 dimers (also referred to herein as "pure OP-1" or "purified OP-1" in contrast to OP-1 as found in media, also referred to herein as soluble or "media" OP-1) have been extensively studied. The conclusion of these 20 studies has been that mature OP-1 is soluble typically only under denaturing conditions. In contrast, the recombinantly produced OP-1 initially secreted into mammalian cell (CHO cell) conditioned media remains soluble in the absence of denaturants. Mature OP-1 has 25 been shown to be soluble in low concentrations of detergents including 0.1% SDS and CHAPS, and in mild denaturing conditions such as low ionic strength at low pH, or in the presence of denaturants with non-ionic detergents, e.g., 6 M urea + 0.3% Tween-80, and in the 30 presence of acidified organic solvent like 50% acetonitrile with 0.1% TFA. It has now been discovered that denaturing solvent conditions separate the pro domain from the mature region, and that previously developed purification protocols involving denaturants

- 18 -

prohibit isolation of the complexed, highly soluble form. As disclosed below, purification of the soluble complexes must be done in the absence of denaturing conditions.

5

10

## 1.b Production of Secreted OP-1 by CHO Cells

Mammalian cell produced OP-1 is synthesized and secreted in a soluble dimeric form, as a result of several post-translational modifications, including appropriate folding, dimerization, glycosylation and cleavage at the juncture of the pro domain and mature domain. Some pro-OP-1 also is secreted without being cleaved, resulting in secreted pro OP-1.

15

#### 1.c Identification of Soluble OP-1 (complexed)

Recombinant OP-1 as expressed by CHO (Chinese

20 hamster ovary cells, see U.S. Patent No. 5,266,638 for
exemplary protocol) is secreted into serum containing
media and exists as in a soluble form. This apparent
solubility of OP-1 could be caused by the association
of OP-1 with a component of the serum or by the

25 secretion of OP-1 from CHO cells in a more soluble form
than its final purified state.

## 1.d Cleveland Mapping of the 39 kDa Protein

30 Using metabolically labelled protein, a methionine-labeled 39 kDa protein from cultured media co-precipitates in an OP-1-dependent manner with OP-1. Protein was further characterized by Cleveland mapping using standard methodologies (see, for example

- 19 -

Cleveland, D.W. (1977) J. Biol. Chem. 252:1102). Bands having an apparent molecular weight of 50, 39 and 19/17 kDa-based on comparison with standard molecular weight markers, were isolated from a PAGE gel, the gel slices 5 placed in the wells of a 20% acrylamide gel along with various amounts of endoproteinase lys-C, electrophoresed into the stacking gel and allowed to digest for 30 min. followed by resolving the generated fragments on the 20% gel. The 50 kDa protein was 10 cleaved to yield two fragments. The larger of these was also generated by the 39 kDa protein while the smaller was generated by the 19/17 kDa proteins. evidence strongly suggested that the 39 kDa protein was the pro domain of OP-1. Further analysis of the secreted form of OP-1 was made possible by the 15 isolation of mature OP-1 as a soluble complex from CHO conditioned media.

# 2. Purification Protocol for Soluble Osteogenic Protein

20

Soluble complexes comprising osteogenic protein can be isolated from conditioned media using a simple, three step chromatographic protocol performed in the absence of denaturants. The protocol involves running the media (or body fluid) over an affinity column, followed by ion exchange and gel filtration chromatographies. The affinity column described below is a Zn-IMAC column. An alternative protocol also envisioned to have utility is an immunoaffinity column, created using standard procedures and, for example, using antibody specific for a given osteogenic protein pro domain (complexed, for example, to a protein Aconjugated Sepharose column.) Protocols for developing

- 20 -

immunoaffinity columns are well described in the art, (see, for example, <u>Guide to Protein Purification</u>, M. Deutscher, ed., Academic Press, San Diego, 1990, particularly sections VII and XI.)

5

In this experiment OP-1 was expressed in CHO cells as described above. The CHO cell conditioned media containing 0.5% FBS was initially purified using Immobilized Metal-Ion Affinity Chromatography (IMAC). The soluble OP-1 complex from conditioned media binds 10 very selectively to the Zn-IMAC resin, presumably through association with the pro domain, and a high concentration of imidazole (50 mM imidazole, pH 8.0) is required for the effective elution of the bound 15 complex. The Zn-IMAC step separates the soluble OP-1 from the bulk of the contaminating serum proteins that elute in the flow through and 35 mM imidazole wash fractions. The Zn-IMAC purified soluble OP-1 is next applied to an S-Sepharose cation-exchange column 20 equilibrated in 20 mM NaPO, (pH 7.0) with 50 mM NaCl. This S-Sepharose step serves to further purify and concentrate the soluble OP-1 complex in preparation for the following gel filtration step. The protein was applied to a Sephacryl S-200HR column equilibrated in Using substantially the same protocol, soluble 25 TBS. osteogenic proteins also may be isolated from one or more body fluids, including serum, cerebro-spinal fluid or peritoneal fluid.

30 IMAC was performed using Chelating-Sepharose (Pharmacia) that had been charged with three column volumes of 0.2 M ZnSO<sub>4</sub>. The conditioned media was titrated to pH 7.0 and applied directly to the ZN-IMAC resin equilibrated in 20 mM HEPES (pH 7.0) with 500 mM

- 21 -

NaCl. The Zn-IMAC resin was loaded with 80 mL of starting conditioned media per mL of resin. After loading the column was washed with equilibration buffer and most of the contaminating proteins were eluted with 35 mM imidazole (pH 7.0) in equilibration buffer. The soluble OP-1 complex is then eluted with 50 mM imidazole (pH 8.0) in 20 mM HEPES and 500 mM NaCl.

The 50 mM imidazole eluate containing the soluble 10 OP-1 complex was diluted with nine volumes of 20 mM  $NaPO_A$  (pH 7.0) and applied to an S-Sepharose (Pharmacia) column equilibrated in 20 mM NaPO (pH 7.0) with 50 mM NaCl. The S-Sepharose resin was loaded with an equivalent of 800 mL of starting conditioned media 15 per mL of resin. After loading the S-Sepharose column was washed with equilibration buffer and eluted with 100 mM NaCl followed by 300 mM and 500 mM NaCl in 20 mM  $NaPO_A$  (pH 7.0). The 300 mM NaCl pool was further purified using gel filtration chromatography. Fifty mls of the 300 mm NaCl eluate was applied to a 5.0 X 90 cm Sephacryl S-200HR (Pharmacia) equilibrated in Tris buffered saline (TBS), 50 mM Tris, 150 mM NaCl (pH 7.4). The column was eluted at a flow rate of 5 mL/minute collecting 10 mL fractions. The apparent 25 molecular of the soluble OP-1 was determined by comparison to protein molecular weight standards (alcohol dehydrogenase (ADH, 150 kDa), bovine serum albumin (BSA, 68 kDa), carbonic anhydrase (CA, 30 kDa) and cytochrome C (cyt C, 12.5 kDa). (see Fig. 3) The purity of the S-200 column fractions was determined by separation on standard 15% polyacrylamide SDS gels stained with coomassie blue. The identity of the mature OP-1 and the pro-domain was determined by N-terminal sequence analysis after separation of the

- 22 -

mature OP-1 from the pro-domain using standard reverse phase C18 HPLC.

Figure 3 shows the absorbance profile at 280 nm.

The soluble OP-1 complex elutes with an apparent molecular weight of 110 kDa. This agrees well with the predicted composition of the soluble OP-1 complex with one mature OP-1 dimer (35-36 kDa) associated with two pro-domains (39 kDa each). Purity of the final complex can be verified by running the appropriate fraction in a reduced 15% polyacrylamide gel.

The complex components can be verified by running the complex-containing fraction from the S-200 or S-15 200HR columns over a reverse phase C18 HPLC column and eluting in an acetonitrile gradient (in 0.1% TFA), using standard procedures. The complex is dissociated by this step, and the pro domain and mature species elute as separate species. These separate species then 20 can be subjected to N-terminal sequencing using standard procedures (see, for example, Guide to Protein Purification, M. Deutscher, ed., Academic Press, San Diego, 1990, particularly pp. 602-613), and the identity of the isolated 36kD, 39kDa proteins 25 confirmed as mature osteogenic protein and isolated, cleaved pro domain, respectively. N-terminal sequencing of the isolated pro domain from mammalian cell produced OP-1 revealed 2 forms of the pro domain, the predominant form being the intact form (beginning 30 at residue 30 of Seq. ID No. 1) and, as a minor species, a truncated form, (beginning at residue 48 of Seg. ID No. 1.) N-terminal sequencing of the polypeptide subunit of the isolated mature species reveals a range of N-termini for the mature sequence,

- 23 -

beginning at residues 293, 300, 313, 315, 316, and 318, of Seq. ID No. 1, all of which are active as demonstrated by the standard bone induction assay. (See, for example, US Pat. Nos. 5,011,691 and 5,266,683 for descriptions of the standard rat bone induction assay.)

## II. Osteogenic Protein Detection

10

As indicated above, the method and compositions of the invention are directed to identifying and/or quantitating preferred forms of osteogenic protein in a solution, such as a culture medium or body fluid. As will be appreciated by those skilled in the art, any means for specifically identifying and quantifying the protein is contemplated. The current state of the art for identifying proteins in solution is by means of an immunoassay, wherein an antibody capable of binding specifically to the protein of interest is used to identify the protein in solution and the amount of bound complex formed then is determined.

Antibody methodologies are well understood and

25 described in the literature. A more detailed
 description of their preparation can be found, for
 example, in Practical Immunology, Butt, W.R., ed.,
 Marcel Dekker, New York, 1984. Broadly, antibodies may
 be raised against one or more preferred forms of an

30 osteogenic protein by immunizing a suitable animal with
 an immunogenic preparation under conditions sufficient
 to induce antibody production in that animal.
 Monoclonal antibodies then can be obtained by fusing
 suitable antibody producing cells such as spleen or

- 24 -

lymph node cells to myeloma cells and screening the fusion products for nuclear reactivity against the immunogen source (e.g., cell line or particular cell type determinant) using standard techniques.

5

The currently preferred method for detecting osteogenic proteins in a solution and/or for quantitating them is by detecting the proteins with osteogenic protein-specific antibodies. The antibodies may be monoclonal or polyclonal in origin, and may be produced by standard methodologies. The osteogenic proteins used as immunogens may be prepared as described above. That is, an intact dimeric species or soluble complex may be used as the antigen (immunogen). Alternatively, a pro domain peptide can be used to advantage, obtained for example, by dissociating a soluble complex and isolating the peptide, or by enzymatic digestion of a precursor form. Alternatively, the entire precursor form may be used as the immunogen.

Antibodies to one or more of these proteins then are raised using standard methods. The antibodies then are exposed to the fluid sample under conditions
25 sufficient to allow specific binding of the antibody to its specific epitope, and the binding partner-osteogenic protein complex formed (here, antibody-osteogenic protein complex) then detected.

30 Immunoassay design considerations include preparation of antibodies (monoclonal or polyclonal) having sufficiently high binding specificity for their antigen that the specifically-bound antibody-antigen complex can be distinguished reliably from nonspecific

interactions. The higher the antibody binding specificity, the lower the antigen concentration that can be detected. The choice of tagging label for detecting osteogenic protein-antibody complex formation 5 also will depend on the detection limitations desired. Enzyme assays (ELISAs) typically allow detection of a colored product formed by interaction of the enzymetagged complex with an enzyme substrate. Alternative labels include radioactive or fluorescent labels. 10 most sensitive label known to date is a chemiluminescent tag where interaction with a reactant results in the production of light. Useful labels include chemiluminescent molecules such as acridium esters or chemiluminescent enzymes where the reactant 15 is an enzyme substrate. When, for example, acridium esters are reacted with an alkaline peroxide solution. an intense flash of light is emitted, allowing the limit of detection to be increased 100 to 10,000 times over those provided by other labels. In addition, the 20 reaction is rapid. A detailed review of chemiluminescence and immunoassays can be found in Weeks, et al., (1983) Methods in Enzymology 133:366-387. Other considerations include the use of microtiter wells or column immunoassays. Column assays 25 may be particularly advantageous where rapidly reacting labels, such as chemiluminescent labels, are used. tagged complex can be eluted to a post-column detector which also contains the reactant or enzyme substrate, allowing the subsequent product formed to be detected 30 immediately.

A detailed review of immunological assay design, theory and protocols can be found in numerous texts in the art, including <u>Practical Immunology</u>, Butt, W.R.,

ed., Marcel Dekker, New York, 1984. Of the various immunoassay formats available, one of the most sensitive is the sandwich technique. In this method. two antibodies capable of binding the analyte of 5 interest are used: one immobilized onto a solid support, and one free in solution, but labeled with some easily detectable chemical compound. As described above, examples of chemical labels that may be used for the second antibody include radioisotopes, fluorescent compounds, and enzymes or other molecules which 10 generate colored or electrochemically active products when exposed to a reactant or enzyme substrate. samples containing analyte (e.g., osteogenic protein in a given form) are placed in this system, the analyte 15 binds to both the immobilized antibody and the labelled antibody. The result is a "sandwich" immune complex on the support's surface. The analyte is detected by washing away nonbound sample components and excess labeled antibody and measuring the amount of labeled 20 antibody complexed to analyte on the support's surface. The sandwich immunoassay is highly specific and very sensitive, provided that labels with good limits of detection are used.

25 Another useful form of immunoassay, particularly useful for screening candidates is the Western blot. Here, proteins of interest are dispersed by gel electrophoresis and immobilized on a nitrocellulose membrane. Candidate antibodies then are added, typically complexed with a means for detection, e.g., radioactive label or enzyme as described above, and complex formation is allowed to occur. The membrane then is washed to remove proteins interacting only by non-specific binding interactions, and complexes that

- 27 -

remain bound are detected. A detailed standard protocol is provided in Molecular Cloning: A Laboratory Manual, Sambrook et al., eds. 2d.ed., Cold Spring Harbor Press, Cold Spring Harbor, 1989. As used berein, typically, osteogenic proteins are electrophoresed under both reducing and oxidizing conditions.

## 2.a Antibody Production

10

Provided below are standard protocols for polyclonal and monoclonal antibody production. For antibodies which recognize the soluble complex form only, preferably the isolated complex itself is used as the antigen. Alternatively, the antigen may comprise the isolated pro domain or a peptide fragment thereof. Where antibodies specific to the mature protein are desired, the antigen preferably comprises the mature dimeric form (e.g., the "purified" form) or a subunit of the dimer comprising at least the C-terminal domain, or a peptide fragment thereof.

# 2a. Polyclonal Antibodies

25 Antibodies then are synthesized as described herein below, and tested <u>in vitro</u> for cross-reactivity with the various forms of the protein of interest.

Polyclonal antibody may be prepared as follows. 30 Each rabbit is given a primary immunization of 100 ug/500  $\mu$ l of antigen, 500  $\mu$ l Complete Freund's Adjuvant. Solubility of a given antigen may be enhanced as needed by combining the antigen in a solubilizing agent, e.g., 0.1% SDS, prior to

- 28 -

combination with the adjuvant. The antigen is injected subcutaneously at multiple sites on the back and flanks of the animal. The rabbit is boosted after a month in the same manner using incomplete Freund's Adjuvant.

5 Test bleeds are taken from the ear vein seven days later. Two additional boosts and test bleeds are performed at monthly intervals until antibody against the osteogenic protein antigen is detected in the serum using an ELISA assay. Then, the rabbit is boosted

10 monthly with 100 µg of antigen and bled (15 ml per bleed) at days seven and ten after boosting.

#### 2b. Monoclonal Antibodies.

Monoclonal antibody specific for a given osteogenic 15 protein may be prepared as follows. A mouse is given two injections of the osteogenic protein antigen. protein or protein fragment preferably is recombinantly produced. The first injection contains 100µg of antigen in complete Freund's adjuvant and is given 20 subcutaneously. The second injection contains 50  $\mu g$  of antigen in incomplete adjuvant and is given intraperitoneally. The mouse then receives a total of 230 ug of OP-1 in four intraperitoneal injections at 25 various times over an extended period of time (e.q., a one to eight month period.) One week prior to fusion, the mouse is boosted intraperitoneally with antiqen (e.g., 100  $\mu$ g) and may be additionally boosted with a pentide fragment conjugated to bovine serum albumin 30 with a suitable cross linking agent. This boost can be repeated five days (IP), four days (IP), three days (IP) and one day (IV) prior to fusion. The mouse spleen cells then are fused to commercially available myeloma cells at a ratio of 1:1 using PEG 1500

- 29 -

(Boehringer Mannheim, Germany), and the fused cells plated and screened for mature or soluble osteogenic protein-specific antibodies using the appropriate portion of the osteogenic protein sequence as antigen.

The cell fusion and monoclonal screening steps readily are performed according to standard procedures well described in standard texts widely available in the art.

## 10 2c. Antibody Specificity.

15

Using these standard procedures, anti-pro domain antisera was prepared from rabbits using the isolated pro domain from OP-1 as the antigen, and monoclonal antibody ("mAb") to the mature domain was produced in mice, using an  $\underline{E.\ coli}$ -produced truncated form of OP-1 as antigen.

Standard Western blot analysis performed as

described herein above under reducing conditions
demonstrates that the anti-pro domain antisera ("antipro") is specific for the pro domain only, while the
mab to mature OP-1 ("anti-mature OP-1") is specific for
the dimer subunits, that the two antibodies do not

cross-react, and that the antibodies and can be used to
distinguish between soluble and mature protein forms in
a sample, e.g., of conditioned media or serum. A
tabular representation of the Western blot results is
in Table I below, where reactivity of mab to mature

OP-1 is indicated by "yy", and reactivity of the
anti-pro antisera is indicated by "xx".

- 30 -

TABLE I

5	Purified Antibody Sol OPl		Conditioned CHO Cell Media	Isolated Pro Domain	Purified Dimer Subunits
	"anti-pro"	xx	xx	xx	
10	"anti- mature OP-1	" уу	уу		уу

In a second series, monoclonal antibodies were

15 raised against each of the following antigens: soluble complex and uncomplexed, mature dimeric species.

Clones then were screened for reactivity against the various forms of OP-1, in an ELISA assay, as

20 described above. Here, the various forms of OP1 tested were immobilized on a surface, the antibody to be screened then was added, and bound antibody detected using a goat anti-mouse antibody. Five different phenotypes or binding categories were identified and

25 are described below. In the table, "S" means soluble complex; "M" means mature, dimeric species, and "P" means isolated pro domain.

- 31 -

#### TABLE II

		м	Protein S	Form P
5	Category	<u> </u>		
	1	+	+	
	2		+	+
	3		+	
10	4			+
	5	+		

Antibodies having the binding character of category #1 recognize an epitope present on both the uncomplexed dimeric species and the soluble form.

Antibodies having the binding character of category #2 recognize an epitope present on both the soluble complex and the pro domain. Antibodies having the binding character of category #3 only recognize the soluble complex form, verifying that a conformational change occurs upon complex formation sufficient to create an epitope not present on other forms of the protein. Antibodies having this binding character are particularly useful for verifying the presence of a complex, including the formation of the soluble complex in vitro from its components (e.g., uncomplexed dimer and isolated pro domain peptide.)

Antibodies having the binding character of category #4 only recognize an epitope on the prodomain, verifying that complex formation is sufficient to mask or destroy an epitope present on the soluble complex. Similarly, antibodies having the binding

- 32 **-**

character of category #5 only recognize an epitope on the mature, dimeric uncomplexed protein form, but not the soluble complex.

Of course, individual members within a given 5 category may bind different epitopes. Accordingly, their binding character with respect to the different protein forms may vary depending on the assay conditions. For example, individual members of category #1, while still recognizing the mature and 10 soluble forms, show preferentially binding affinity for the soluble form over the mature form under sandwich ELISA conditions where the member of category #1 constituted the capture antibody. Similarly, individual members of category #2 demonstrate variable 15 binding for the pro domain under Western blot conditions.

#### III. <u>Immunoassays</u>

20

The ability to detect osteogenic proteins in solution and to distinguish between soluble and mature dimeric forms provides a valuable tool for protein production systems. Quality control considerations

25 require that means be available for determining both the form of the protein in solution and its quantity. This is particularly true for biological therapeutics where pharmacologically certain forms must be provided for clinical use. The method also provides a useful tool for diagnostic assays, allowing one to monitor the level and type of osteogenic protein free in the body, e.g., in serum and other body fluids.

- 33 -

A currently preferred detection means for evaluating the level of osteogenic protein in a fluid, including culture or a body fluid, comprises an immunoassay utilizing an antibody or other suitable 5 binding protein capable of reacting specifically with an osteogenic protein and being detected as part of a complex with the protein. Immunoassays may be performed using standard techniques known in the art and using antibodies raised against the protein and 10 specific for that protein.

Antibodies which recognize an osteogenic protein form of interest may be generated as described herein and these antibodies then used to monitor the levels of protein in a fluid, including a body fluid, such as serum, whole blood or peritoneal fluid.

To monitor endogenous concentrations of the soluble form of the protein, the antibody chosen preferably has 20 binding specificity for the soluble form. For endogenous proteins, these antibodies may have specificity for the pro domain and/or the soluble complex (e.g., binding categories 1-3, above). Such antibodies may be generated by using the pro domain or 25 a portion thereof as the antigen, or the soluble complex itself, essentially as described herein. suitable pro domain for use as an antigen may be obtained by isolating the soluble complex and then separating the noncovalently associated pro domain from 30 the mature domain using standard procedures, e.g., by chromatographic means (preferably using an ion-exchange column under denaturing conditions, e.g., 6M urea), as described above or by separation by gel electrophoresis. Alternatively, the pro form of the

- 34 -

protein in its monomeric form may be used as the antigen and the candidate antibodies screened by Western blot or other standard immunoassay for those which recognize the pro domain of the soluble form of the protein of interest, but not the mature form, also as described above.

Monomeric pro forms can be obtained from cell lysates of CHO produced cells, or from prokaryotic expression of a DNA encoding the pro form, in for example, <a href="E.coli">E.coli</a>, or from a commercially available bacculovirus expression system in insect cells. The pro form, which has an apparent molecular weight of about 50 kDa in mammalian cells, can then be isolated as described above.

In order to detect and/or quantitate the amount of osteogenic protein present in a solution, an immunoassay may be performed to detect the osteogenic 20 protein using a polyclonal or monoclonal antibody specific for that protein. Here, soluble and mature forms of the osteogenic protein also may be distinguished by using antibodies that discriminate between the two forms of the proteins as described 25 above. Currently preferred assays include ELISAS and radioimmunassays, including standard competitor assays useful for quantitating the osteogenic protein in a sample, where an unknown amount of sample protein is allowed to react with anti-osteogenic protein antibody and this interaction is competed with a known amount of 30 labeled antigen. The level of bound or free labeled antigen at equilibrium then is measured to quantitate the amount of unlabeled antigen in solution, the amount of sample antigen being proportional to the amount of

- 35 -

free labeled antigen. Exemplary protocols for these assays are provided below. However, as will be appreciated by those skilled in the art, variations of these protocols, as well as other immunoassays, 5 including Western blots, are well known in the literature and within the skill of the art. For example, in the ELISA protocol provided below, soluble OP-1 is identified in a sample using biotinylated antipro antiserum. Biotinylated antibodies can be 10 visualized in a colormetric assay or in a chemiluminescent assay, as described below. Alternatively, the antibody can be radio-labeled with a suitable molecule, such as 125 I. Still another protocol that may be used is a solid phase immunoassay. 15 preferably using an affinity column with antiosteogenic protein antibody complexed to the matrix surface and over which a serum sample may be passed. detailed description of useful immunoassays, including protocols and general considerations is provided in, 20 for example, Molecular Cloning: A Laboratory Manual, Sambrook et al., eds. Cold Spring Harbor Press, New York, 1989, particularly Section 18.

For serum assays, the serum preferably first is
partially purified to remove some of the excess,
contaminating serum proteins, such as serum albumin.
Preferably the serum is extracted by precipitation in
ammonium sulfate (e.g., 45%) such that the complex is
precipitated. Further purification can be achieved
using purification strategies that take advantage of
the differential solubility of soluble osteogenic
protein complex or mature osteogenic proteins relative
to that of the other proteins present in serum.

- 36 -

Further purification also can be achieved by chromatographic techniques well known in the art.

## 3a. Assays

5

Soluble OP-1 may be detected using a polyclonal or monoclonal antibody in an ELISA, as described below in this experiment, polyclonal antibody specific for the OP-1 pro domain is omalyzed1 ug/100 ul of 10 affinity-purified polyclonal rabbit IgG specific for OP-1-pro is added to each well of a 96-well plate and incubated at 37°C for an hour. The wells are washed four times with 0.167M sodium borate buffer with 0.15 M NaCl (BSB), pH 8.2, containing 0.1% Tween 20. 15 minimize non-specific binding, the wells are blocked by filling completely with 1% bovine serum albumin (BSA) in BSB and incubating for 1 hour at 37°C. are then washed four times with BSB containing 0.1% Tween 20. A 100  $\mu$ l aliquot of an appropriate dilution 20 of each of the test samples of cell culture supernatant or serum sample is added to each well in triplicate and incubated at 37°C for 30 min. After incubation, 100 µl biotinylated antibody from rabbit anti-pro anti-serum (stock solution is about 1 mg/ml and diluted 1:400 in 25 BSB containing 1% BSA before use) is added to each well and incubated at 37°C for 30 min. The wells are then washed four times with BSB containing 0.1% Tween 20. 100 µl strepavidin-alkaline (Southern Biotechnology Associates, Inc. Birmingham, Alabama, diluted 1:2000 in 30 BSB containing 0.1% Tween 20 before use) is added to each well and incubated at 37°C for 30 min. are washed four times with 0.5M Tris buffered Saline (TBS), pH 7.2. 50µl substrate (ELISA Amplification System Kit, Life Technologies, Inc., Bethesda, MD) is

- 37 -

added to each well incubated at room temperature for 15 Then, 50  $\mu$ l amplifier (from the same amplification system kit) is added and incubated for another 15 min at room temperature. The reaction is 5 stopped by the addition of 50  $\mu$ l 0.3 M sulphuric acid. The OD at 490 nm of the solution in each well is recorded. To quantitate the level of soluble OP-1 in the sample, a standard curve is performed in parallel with the test samples. In the standard curve, known 10 increasing amounts of purified OP-1-pro is added. Alternatively, using, for example, Lumi-phos 530 (Analytical Luminescence Laboratories) as the substrate and detection at 300-650 nm in a standard luminometer, complexes can be detected by chemiluminescence, which 15 typically provides a more sensitive assay than detection by means of a visible color change.

## 3b. Plate-based Radioimmunoassay

20 Osteogenic protein (soluble or mature form) may be detected in a standard plated-based radioimmunoassay as follows. Empirically determined limiting levels of anti-osteogenic protein antibody (e.g., anti-OP-1, typically 50-80 ng/well) are bound to wells of a PVC 25 plate e.g., in 50 µl PBS phosphate buffered saline. After sufficient incubation to allow binding at room temperature, typically one hour, the plate is washed in a borate-buffered saline/Tween 20 solution, ("washing buffer"), and 200  $\mu$ l of block (3% BSA, 0.1M lysine in 30 1xBSB) is added to each well and allowed to incubate for 1 hour, after which the wells are washed again in washing buffer. 40  $\mu$ l of a sample composed of serially diluted plasma (preferably partially purified as described above) or osteogenic protein standard (e.g.,

- 38 -

OP-1) is added to wells in triplicate. Samples preferably are diluted in PTTH (15 mM KH<sub>2</sub>PO<sub>4</sub>, 8 mM Na<sub>2</sub>PO<sub>4</sub>, 27 mM KCl, 137 mM NaCl, 0.05% Tween 20, 1 mg/ml HSA, 0.05% NaN<sub>3</sub>, pH 7.2). 10 µl of labelled competitor antigen, preferably 100,000-500,000 cpm/sample is added (e.g., 125 OP-1, radiolabelled using standard procedures), and plates are incubated overnight at 4°C. Plates then are washed in washing buffer, and allowed to dry. Wells are cut apart and bound labelled OP-1 counted in a standard gamma counter. The quantities of bound labelled antigen (e.g., 125 OP-1) measured in the presence and absence of sample then are compared, the difference being proportional to the amount of sample antigen (osteogenic protein) present in the

## 3c. Production Monitoring Considerations

Samples for testing the level of protein production
includes culture supernatants or cell lysates,
collected periodically and evaluated for OP-1
production by immunoblot analysis (Sambrook et al.,
eds., 1989, Molecular Cloning, Cold Spring Harbor
Press, Cold Spring Harbor, NY), or a portion of the
cell culture itself, collected periodically and used to
prepare polyA+ RNA for mRNA analysis. To monitor de
novo OP-1 synthesis, some cultures are labeled
according to conventional procedures with an

35 S-methionine/35 S-cysteine mixture for 6-24 hours and
then evaluated to OP-1 synthesis by conventional
immunoprecipitation methods.

- 39 -

# 3.d <u>Diagnostics using Antibodies to Soluble</u> Osteogenic <u>Protein Complex</u>

The antibodies of this invention also may be used to monitor the level of soluble protein in the body. Fluctuations in osteogenic protein levels present in the bloodstream or peritoneal fluid then may be used to evaluate tissue viability. For example, osteogenic proteins are detected associated with regenerating tissue and/or may be released from dying cells into surrounding peritoneal fluid.

Serum samples may be obtained by standard venipuncture and serum prepared by centrifugation at 3,000 RPM for ten minutes. Similarly, peritoneal fluid samples may be obtained by a standard fluid extraction methodology. The presence of osteogenic protein in the serum or peritoneal fluid then may be assessed by standard Western blot (immunoblot), ELISA or RIA 20 procedures. Briefly, for example, with the ELISA, samples may be diluted in an appropriate buffer, such as phosphate-buffered saline, and 50 µl aliquots allowed to absorb to flat bottomed wells in microtitre plates pre-coated with soluble osteogenic protein-25 specific antibody, and allowed to incubate for 18 hours at 4°C. Plates then may be washed with a standard buffer and incubated with 50 µl aliquots of a second osteogenic protein-specific antibody conjugated with a detecting agent, e.g., biotin, in an appropriate buffer, for 90 minutes at room temperature. Osteogenic 30 protein-antibody complexes then may be detected using standard procedures.

- 40 -

Alternatively, an osteogenic protein-specific affinity column may be created using, for example, soluble osteogenic protein-specific antibodies adsorbed to a column matrix, and passing the fluid sample 5 through the matrix to selectively extract the protein of interest. The protein then is eluted. A suitable elution buffer may be determined empirically by determining appropriate binding and elution conditions first with a control (e.g., purified, 10 recombinantly-produced protein.) Fractions then are tested for the presence of the soluble form protein by standard immunoblot. Protein concentrations in serum or other fluid samples then may be determined using standard protein quantification techniques, including 15 by spectrophotometric absorbance or by quantitation by ELISA or RIA antibody assays. Using this procedure, OP-1 has been identified in serum.

OP-1 was detected in human serum using the 20 following assay. A monoclonal antibody raised against mammalian, recombinantly produced OP-1 using standard immunology techniques well described in the art and described generally herein, was immobilized by passing the antibody over an activated agarose gel (e.g., Affi-25 Gel<sup>TM</sup>, from Bio-Rad Laboratories, Richmond, CA, prepared following manufacturer's instructions), and used to purify OP-1 from serum. Human serum then was passed over the column and eluted with 3M K-thiocyanate. K-thiocyanante fractions then were 30 dialyzed in 6M urea, 20mM PO, pH 7.0, applied to a C8 HPLC column, and eluted with a 20 minute, 25-50% acetonitrile/0.1% TFA gradient. Mature, recombinantly produced OP-1 homodimers elute between 20-22 minutes. Accordingly, these fractions from the affinity-purified

- 41 -

human serum sample were collected and tested for the presence of OP-1 by standard immunoblot using an OP-1-specifc antibody, and the protein identity confirmed by N-terminal sequencing.

5

The invention may be embodied in other specific forms without departing from the spirit or essential characteristics thereof. The present embodiments are therefore to be considered in all respects as illustrative and not restrictive, the scope of the invention being indicated by the appended claims rather than by the foregoing description, and all changes which come within the meaning and range of equivalency of the claims are therefore intended to be embraced

15 therein.

- 42 -

## SEQUENCE LISTING

5	(1)	GENERAL	INFORMATION:
10 15		(i)	APPLICANT: (A) NAME: CREATIVE BIOHOLECULES, INC. (B) STREET: 45 SOUTH STREET (C) CITY: HOPKINTON (D) STATE: HA (E) COUNTRY: USA (F) POSTAL CODE (ZIP): 01748 (G) TELEPHONE: 508-435-9001 (H) TELEFAX: 508-435-0992 (I) TELEX:
		(ii)	TITLE OF INVENTION: METHODS AND COMPOSITIONS FOR RECOMBINANT OSTEOGENIC PROTEIN PRODUCTION
20		(iii)	NUMBER OF SEQUENCES: 6
25		(iv)	CORRESPONDENCE ADDRESS:  (A) ADDRESSEE: TESTA, HURWITZ & THIBEAULT  (B) STREET: 53 STATE STREET  (C) CITY: BOSTON  (D) STATE: HASSACHUSETTS  (E) COUNTRY: U.S.A.  (F) ZIP: 02109
30		(▽)	COMPUTER READABLE FORM:  (A) MEDIUM TYPE: Floppy disk  (B) COMPUTER: IBM PC compatible  (C) OPERATING SYSTEM: PC-DOS/MS-DOS  (D) SOFTWARE: Patent In Release #1.0, Version #1.25
35			(b) bollward. Ideate in McLease with velocin with
		(vi)	CURRENT APPLICATION DATA:  (A) APPLICATION NUMBER: US  (B) FILING DATE: HEREWITH  (C) CLASSIFICATION:
40		(vii)	ATTORNEY/AGENT INFORMATION: (A) NAME: PITCHER, EDMUND R. (B) REGISTRATION NUMBER: 27,829 (C) REFERENCE/DOCKET NUMBER: CRP-096PC
45		(viii)	TELECOMMUNICATION INFORMATION: (A) TELEPHONE: 617/248-7000 (B) TELEFAX: 617/248-7100

5	(2)	IN	FORM.	) SI (4 (1	EQUENA) 1 B) 1 C) 1	R SEC NCE ( LENGT LYPE: STRAI LOPOI	CHARA  CH:  nuc  NDEDI	ACTE 1822 cleic NESS	RIST: base ac: si:	e pa: id	irs						
10			(ii	) <b>M</b> (	DLEC	ULE :	TYPE:	: cDl	AN								
15			(vi	(1 (1	A) ( F) :	NAL S ORGAI	NISM:	: HO									
20			(ix	( <i>I</i>	B) 1	NAME /	PION: R IN: /pro/ ev:	6 49 FORM Oducti Ldeno	134 ATION t= "N		-PP" RIMEN	NTAL	= "0	STEO	GENI(	C PROTEIN"	
25	GGT	CGG(	(xi)	·	-	NCE I				•					t His	C GTG s Val	57
30	CGC Arg	TCA Ser 5	CTG Leu	CGA Arg	GCT Ala	GCG Ala	GCG Ala 10	CCG Pro	CAC His	AGC Ser	TTC Phe	GTG Val 15	GCG Ala	CTC Leu	TGG Trp	GCA Ala	105
35														CTG Leu			153
40														CAG Gln			201
45														CCC Pro 65			249
40														ATG Het			297
50	CTC	CAC	CTC	ТΔС	AAC	CCC	ATC	ccc	стс	GAG	CAC	ccc	ccc	ccc	ccc	ccc	345

- 44 -

	Leu	Asp 85	Leu	Tyr	Asn	Ala	Met 90	Ala	Val	Glu	Glu	Gly 95	Gly	Gly	Pro	Gly	
5	GGC Gly 100	CAG Gln	GGC Gly	TTC Phe	TCC Ser	TAC Tyr 105	CCC Pro	TAC Tyr	AAG Lys	GCC Ala	GTC Val 110	TTC Phe	AGT Ser	ACC Thr	CAG Gln	GGC Gly 115	393
10												CTC Leu					441
			Met									GAC Asp					489
15												GAT Asp					537
20												CGG Arg 175					585
25												CGG Arg					633
30	CAG Gln	GTG Val	CTC Leu	CAG Gln	GAG Glu 200	CAC His	TTG Leu	GGC Gly	AGG Arg	GAA Glu 205	TCG Ser	GAT Asp	CTC Leu	TTC Phe	CTG Leu 210	CTC Leu	681
			Arg									TGG Trp					729
35												CCG Pro					777
40	GGC Gly	CTG Leu 245	CAG Gln	CTC Leu	TCG Ser	GTG Val	GAG Glu 250	ACG Thr	CTG Leu	GAT Asp	GGG Gly	CAG Gln 255	AGC Ser	ATC Ile	AAC Asn	CCC Pro	825
45												CAG Gln					873
50												CAC His					921

															ACG Thr			969
5	AAG Lys														AGC Ser		1	1017
10															AGC Ser		1	1065
15	CGA Arg 340	GAC Asp	CTG Leu	GGC Gly	TGG Trp	CAG Gln 345	GAC Asp	TGG Trp	ATC Ile	ATC Ile	GCG Ala 350	CCT Pro	GAA Glu	GGC Gly	TAC Tyr	GCC Ala 355	1	.113
20	GCC Ala	TAC Tyr	TAC Tyr	TGT Cys	GAG Glu 360	GGG Gly	GAG Glu	TGT Cys	GCC Ala	TTC Phe 365	CCT Pro	CTG Leu	AAC Asn	TCC Ser	TAC Tyr 370	ATG Met	1	.161
20	AAC Asn	GCC Ala	ACC Thr	AAC Asn 375	CAC His	GCC Ala	ATC Ile	GTG Val	CAG Gln 380	ACG Thr	CTG Leu	GTC Val	CAC His	TTC Phe 385	ATC Ile	AAC Asn	1	.209
25	CCG Pro	GAA Glu	ACG Thr 390	GTG Val	CCC Pro	AAG Lys	CCC Pro	TGC Cys 395	TGT Cys	GCG Ala	CCC Pro	ACG Thr	CAG Gln 400	CTC Leu	AAT Asn	GCC Ala	1	.257
30	ATC Ile	TCC Ser 405	GTC Val	CTC Leu	TAC Tyr	TTC Phe	GAT Asp 410	GAC Asp	AGC Ser	TCC Ser	AAC Asn	GTC Val 415	ATC Ile	CTG Leu	AAG Lys	AAA Lys	1	.305
35				ATG Met									TAGO	TCCI	rcc		1	351
	GAGA	ATTO	CAG A	CCCI	TTGG	G GC	CAAC	TTTT	TCI	GGAT	CCT	CCAT	TGCI	CG C	CTTG	GCCAG	. 1	411
40	GAAC	CAGO	AG A	CCAA	CTGC	C TI	TTG1	GAGA	CCI	TCCC	стс	CCTA	TCCC	CA A	CTTI	'AAAGG	. 1	471
40	TGTG	AGAG	TA I	TAGG	AAAC	A TO	AGCA	GCAT	' ATG	GCTI	TTG	ATCA	GTTI	TT C	CAGTG	GCAGC	1	531
	ATCC	CTAA	AA C	CAAGA	TCCI	'A CA	AGCI	GTGC	AGG	CAAA	ACC	TAGO	AGGA	AA A	AAAA	ACAAC	1	591
45	GCAT	'AAAC	AA A	AAATG	GCCG	G GC	CAGG	TCAT	TGG	CTGG	GAA	GTCI	CAGC	CA T	CGCAC	GGACT	1	651
	CGTT	TCCA	GA C	GTAA	CATTA	G AC	CGCC	TACC	AGC	CAGG	CCA	CCCA	GCCG	TG C	GAGG	AAGGG	1	711
50	GGCG	TGGC	CAA C	GGGT	rggg	CA CA	TTGO	TGT	TG1	GCGA	AAG	GAAA	ATTO	AC (	CCGGA	AGTTC	1	771
	CTGT	CAATA	L AA	CGTCA	CAAT	CA AA	ACGA	ATGA	ATO	AAAA	AAA	AAAA	AAAA	AA A	A.		1	822

PCT/US94/02335

- 46 -

	(2)	TNI	UKM	YIIOD	ı rur	SEC	ί τη	NO: 2	• •							
5			(i)	) SE (# (F	i) I	ENGT	CHARA CH: 4 : ami LOGY:	31 a	umino acid		ids					
10			<b>(ii</b> )	) MC	LECU	JLE 1	TYPE:	pro	teir	1						
			(ix)		EATUE O) O		INFO	ORMAT	CION:	: /1	?rodı	ıct='	'OP1-	-PP"		
15			(xi	) SI	EQUE	ICE I	DESCI	RIPT	ON:	SEQ	ID 1	10:2:				
	Met 1	His	Val	Arg	Ser 5	Leu	Arg	Ala	Ala	Ala 10	Pro	His	Ser	Phe	Val 15	Ala
20	Leu	Trp	Ala	Pro 20	Leu	Phe	Leu	Leu	Arg 25	Ser	Ala	Leu	Ala	Asp 30	Phe	Ser
25	Leu	Asp	Asn 35	Glu	Val	His	Ser	Ser 40	Phe	Ile	His	Arg	Arg 45	Leu	Arg	Ser
25	Gln	Glu 50	Arg	Arg	Glu	Met	Gln 55	Arg	Glu	Ile	Leu	Ser 60	Ile	Leu	Gly	Leu
30	Pro 65	His	Arg	Pro	Arg	Pro 70	His	Leu	Gln	Gly	Lys 75	His	Asn	Ser	Ala	Pro 80
	Met	Phe	Het	Leu	Asp 85	Leu	Tyr	Asn	Ala	Met 90	Ala	Val	Glu	Glu	Gly 95	Gly
35	Gly	Pro		Gly 100	<b>Gl</b> n	Gly	Phe	Ser	Tyr 105	Pro	Tyr	Lys	Ala	Val 110	Phe	Ser
	Thr	Gln	Gly 115	Pro	Pro	Leu	Ala	Ser 120	Leu	G1n	Asp	Ser	His 125	Phe	Leu	Thr
40	Asp	Ala 130	Asp	Met	Val	Met	Ser 135	Phe	Val	Asn	Leu	Val 140	Glu	His	Asp	Lys
45	Glu 145	Phe	Phe	His	Pro	Arg 150	Tyr	His	His	Arg	Glu 155	Phe	Arg	Phe	Asp	Leu 160
	Ser	Lys	Ile	Pro	Glu 165	Gly	Glu	Ala	Val	Thr 170	Ala	Ala	Glu	Phe	Arg 175	Ile
50	Tyr	Lys		Tyr 180	Ile	Arg	Glu	Arg	Phe 185	Asp	Asn	Glu	Thr	Phe 190	Arg	Ile

- 47 -

	Ser	Val	<b>Tyr</b> 195	Gln	Val	Leu	Gln	Glu 200	His	Leu	Gly	Arg	G1u 205	Ser	Asp	Leu
5	Phe	Leu 210	Leu	Asp	Ser	Arg	Thr 215	Leu	Trp	Ala	Ser	Glu 220	Glu	Gly	Trp	Leu
10	Val 225	Phe	Asp	Ile	Thr	Ala 230	Thr	Ser	Asn	His	Trp 235	Val	Val	Asn	Pro	Arg 240
10	His	Asn	Leu	Gly	Leu 245	Gln	Leu	Ser	Val	Glu 250	Thr	Leu	Asp	Gly	Gln 255	Ser
15	Ile	Asn	Pro	Lys 260	Leu	Ala	Gly	Leu	Ile 265	Gly	Arg	His	Gly	Pro 270	Gln	Asn
	Lys	Gln	Pro 275	Phe	Het	Val	Ala	Phe 280	Phe	Lys	Ala	Thr	Glu 285	Val	His	Phe
20	Arg	Ser 290	Ile	Arg	Ser	Thr	Gly 295	Ser	Lys	Gln	Arg	Ser 300	Gln	Asn	Arg	Ser
25	Lys 305	Thr	Pro	Lys	Asn	Gln 310	Glu	Ala	Leu	Arg	Met 315	Ala	Asn	Val	Ala	Glu 320
	Asn	Ser	Ser	Ser	Asp 325	Gln	Arg	Gln	Ala	Cys 330	Lys	Lys	His	Glu	Leu 335	Tyr
30	Val	Ser		Arg 340	Asp	Leu	Gly	Trp	Gln 345	Asp	Trp	Ile	Ile	Ala 350	Pro	Glu
	Gly	Tyr	Ala 355	Ala	Tyr	Tyr	Cys	Glu 360	Gly	Glu	Cys	Ala	Phe 365	Pro	Leu	Asn
35	Ser	Tyr 370	Met	Asn	Ala	Thr	Asn 375	His	Ala	Ile	Val	Gln 380	Thr	Leu	Val	His
40	Phe 385	Ile	Asn	Pro	Glu	Thr 390	Val	Pro	Lys	Pro	Cys 395	Cys	Ala	Pro	Thr	Gln 400
40	Leu	Asn	Ala	Ile	Ser 405	Val	Leu	Tyr	Phe	Asp 410	Asp	Ser	Ser	Asn	Val 415	Ile
45	Leu	Lys		Tyr 420	Arg	Asn	Met	Val	Val 425	Arg	Ala	Cys	Gly	Cys 430	His	
	(2)	INF	ORMA:	CION	FOR	SEQ	ID 1	10:3:							•	
50		(i)	(4	À) LI	ENGTI	iARAC i: 17 nucl	723 l	oase	pair	rs						

- 48 -

(C) STRANDEDNESS: single (D) TOPOLOGY: linear

_	(11) MOLECULE TYPE: cDNA	
5	(vi) ORIGINAL SOURCE: (A) ORGANISM: Homo sapiens (F) TISSUE TYPE: HIPPOCAMPUS	
10	(ix) FEATURE:  (A) NAME/KEY: CDS  (B) LOCATION: 4901696  (D) OTHER INFORMATION: /function= "OSTEOGENIC PROTEIN"	
15	/product= "hOP2-PP" /note= "hOP2 (cDNA)"	
	(xi) SEQUENCE DESCRIPTION: SEQ ID NO:3:	
20	GGCGCCGGCA GAGCAGGAGT GGCTGGAGGA GCTGTGGTTG GAGCAGGAGG TGGCACGGCA	60
	GGGCTGGAGG GCTCCCTATG AGTGGCGGAG ACGGCCCAGG AGGCGCTGGA GCAACAGCTC	120
	CCACACCGCA CCAAGCGGTG GCTGCAGGAG CTCGCCCATC GCCCCTGCGC TGCTCGGACC	180
25	GCGGCCACAG CCGGACTGGC GGGTACGGCG GCGACAGAGG CATTGGCCGA GAGTCCCAGT	240
	CCGCAGAGTA GCCCCGGCCT CGAGGCGGTG GCGTCCCGGT CCTCTCCGTC CAGGAGCCAG	300
30	GACAGGTGTC GCGCGGCGG GCTCCAGGGA CCGCGCCTGA GGCCGGCTGC CCGCCCGTCC	360
	CGCCCCGCCC CGCCGCCGC CGCCCGCCGA GCCCAGCCTC CTTGCCGTCG GGGCGTCCCC	420
	AGGCCCTGGG TCGGCCGCGG AGCCGATGCG CGCCCGCTGA GCGCCCCAGC TGAGCGCCCC	480
35	CGGCCTGCC ATG ACC GCG CTC CCC GGC CCG CTC TGG CTC CTG GGC CTG  Met Thr Ala Leu Pro Gly Pro Leu Trp Leu Leu Gly Leu  1 5 10	528
40	GCG CTA TGC GCG CTG GGC GGG GGC CCC GGC CTG CGA CCC CCG CCC Ala Leu Cys Ala Leu Gly Gly Gly Gly Pro Gly Leu Arg Pro Pro Pro 15 20 25	576
45	GGC TGT CCC CAG CGA CGT CTG GGC GCG CGC GAG CGC CGG GAC GTG CAG Gly Cys Pro Gln Arg Arg Leu Gly Ala Arg Glu Arg Arg Asp Val Gln 30 40 45	624
50	CGC GAG ATC CTG GCG GTG CTC GGG CTG CCT GGG CGG C	672

- 49 -

		CCC Pro							720
5		CTG Leu 80							768
10		GAG Glu							816
15		GTG Val							864
20		TTC Phe							912
		GCG Ala							960
25		ACC Thr 160							1008
30		GAG Glu							1056
35		GAG Glu							1104
40		CTG Leu							1152
		GAC Asp							1200
45		GCC Ala 240							1248
50		CCG Pro							1296

- 50 -

5		AGG Arg															1344
,		GGG Gly															1392
10		CGG Arg															1440
<b>1</b> 5	TGG Trp	GTC Val	ATC Ile 320	GCT Ala	CCC Pro	CAA Gln	GGC Gly	TAC Tyr 325	TCG Ser	GCC Ala	TAT Tyr	TAC Tyr	TGT Cys 330	GAG Glu	GGG Gly	GAG Glu	1488
20	TGC Cys	TCC Ser 335	TTC Phe	CCA Pro	CTG Leu	GAC Asp	TCC Ser 340	TGC Cys	ATG Met	AAT Asn	GCC Ala	ACC Thr 345	AAC Asn	CAC His	GCC Ala	ATC Ile	1536
25	CTG Leu 350	CAG Gln	TCC Ser	CTG Leu	GTG Val	CAC His 355	CTG Leu	ATG Met	AAG Lys	CCA Pro	AAC Asn 360	GCA Ala	GTC Val	CCC Pro	AAG Lys	GCG Ala 365	1584
23		TGT Cys															1632
30	AGC Ser	AGC Ser	AAC Asn	AAC Asn 385	GTC Val	ATC Ile	CTG Leu	CGC Arg	AAA Lys 390	GCC Ala	CGC Arg	AAC Asn	ATG Met	GTG Val 395	GTC Val	AAG Lys	1680
35		TGC Cys				T GA	AGTC	AGCCC	GCC	CCAGO	CCT	ACTO	CAG				1723
40	(2)	INFO		TION SEQUE													
45		`	, -	(A) (B)	LEN TYI	IGTH:	402 umino GY: ]	ami aci	ino a id		3						
4.7				IOLEC EQUE						) ID	NO: 4	:			•		
50	Met 1	Thr	Ala	Leu	Pro 5	Gly	Pro	Leu	Trp	Leu 10	Leu	Gly	Leu	Ala	Leu 15	Cys	

	AIA	Leu	GIÀ	20	GIŻ	Gly	Pro	GIY	Leu 25	Arg	Pro	Рто	Pro	30	Cys	Pro
5	Gln	Arg	Arg 35	Leu	Gly	Ala	Arg	Glu 40	Arg	Arg	Asp	Val	Gln 45	Arg	Glu	Ile
10	Leu	Ala 50	Val	Leu	Gly	Leu	Pro 55	Gly	Arg	Pro	Arg	Pro 60	Arg	Ala	Pro	Pro
	Ala 65	Ala	Ser	Arg	Leu	Pro 70	Ala	Ser	Ala	Pro	Leu 75	Phe	Het	Leu	Asp	Leu 80
15	Tyr	His	Ala	Met	Ala 85	Gly	Asp	Asp	Asp	Glu 90	Asp	Gly	Ala	Pro	Ala 95	Glu
	Arg	Arg	Leu	Gly 100	Arg	Ala	Asp	Leu	Val 105	Ket	Ser	Phe	Val	Asn 110	Met	Val
20	Glu	Arg	Asp 115	Arg	Ala	Leu	Gly	His 120	Gln	Glu	Pro	His	Trp 125	Lys	Glu	Phe
25	Arg	Phe 130	Asp	Leu	Thr	Gln	Ile 135	Pro	Ala	Gly	Glu	Ala 140	Val	Thr	Ala	Ala
	Glu 145	Phe	Arg	Ile	Tyr	Lys 150	Val	Pro	Ser	Ile	His 155	Leu	Leu	Asn	Arg	Thr 160
30					165					170				Asn	175	
				180					185					Gly 190	_	
35	Gly	Trp	Leu 195	Val	Leu	Asp	Val	Thr 200	Ala	Ala	Ser	Asp	Cys 205	Trp	Leu	Leu
40		210					215				-	220		Thr		•
	225					230					235			Gln		240
45					245					250			_	Ala	255	
				260			•		265				-	Arg 270	_	
50	Pro	Lys	Lys 275	Ser	Asn	Glu	Leu	Pro 280	Gln	Ala	Asn	Arg	Leu 285	Pro	Gly	Ile

- 52 -

	Phe	Asp 290	Asp	Val	His	Gly	Ser 295	His	Gly	Arg	Gln	Val 300	Cys	Arg	Arg	His	
5	Glu 305	Leu	Tyr	Val	Ser	Phe 310	Gln	Asp	Leu	Gly	Trp 315	Leu	Asp	Trp	Val	11e 320	
10	Ala	Pro	Gln	Gly	Tyr 325	Ser	Ala	Tyr	Tyr	Cys 330	Glu	Gly	Glu	Cys	Ser 335	Phe	
10	Pro	Leu	Asp	Ser 340	Cys	Met	Asn	Ala	Thr 345	Asn	His	Ala	Ile	Leu 350	Gln	Ser	
15	Leu	Val	His 355	Leu	Met	Lys	Pro	Asn 360	Ala	Val	Pro	Lys	Ala .365	Cys	Cys	Ala	
	Pro	Thr 370	Lys	Leu	Ser	Ala	Thr 375	Ser	Val	Leu	Туг	Туг 380	Asp	Ser	Ser	Asn	
20	Asn 385	Val	Ile	Leu	Arg	Lys 390	Ala	Arg	Asn	Met	Val 395	Val	Lys	Ala	Cys	Gly 400	
	Cys	His															
25	(2)	INF	ORMA:	rion	FOR	SEQ	ID I	NO:5	:								
30		<b>(i</b> )	(4 (1	QUENC A) LI B) T D) T	ENGTI YPE:	H: 10 amin	02 au	mino cid		ds							
		<b>(ii</b> )	) но	LECU	LE T	YPE:	pro	tein									
35		(ix	. (A	ATURI A) NA B) L(	AME/I												
40			(I	o) o:	/ne FR	ote= OM A	"WH	EREII UP O	N EA	CH X. E OR	AA I MOR	S IN	ECIF			SELEC	
45		(xi	) SE	QUEN	CE D	ESCR:	IPTI	ON:	SEQ :	ID N	0:5:						
45		Cy:	s Xa	a Xa	a Hi	5 G1:	ı Le	и <b>Ту</b> :	r Va	l Xa	a Ph 10		a As	p Le	u Gly	y Trp 15	Xaa
50		As	p Tr	р Ха	a Il 20	e Al	a Pr	о Ха	a Gl	y Ty 25		a Al	а Ту	т Ту	r Cy 30	s Glu	Gly

	Glu	Cys	<b>Xaa</b> 35	Phe	Pro	Leu	Xaa	Ser 40	Xaa	Met	Asn	Ala	Thr 45	Asn	His	Ala
5	Ile	Xaa 50	Gln	Xaa	Leu	Val	His 55	Xaa	Xaa	Xaa	Pro	Xaa 60	Xaa	Val	Pro	Lys
	Xaa 65	Cys	Cys	Ala	Pro	Thr 70	Xaa	Leu	Xaa	Ala	Xaa 75	Ser	Val	Leu	Tyr	Xaa 80
10	Asp	Xaa	Ser	Xaa	Asn 85	Val	Xaa	Leu	Xaa	Lys 90	Xaa	Arg	Asn	Met	Val 95	Val
15	Xaa	Ala	Cys	Gly 100	Cys	His										
13	(ix)	(A)	) NAI	: ME/KI CATI(							•					
20				HER : /not FROI	INFOI te= ' I A (	RMAT: WHEI ROUL	CON: REIN POF	EACH ONE	OR I	OPX A IS MORE FICAT	SPE	CIFIE	ED AL	ONI	ACII	os
25	(xi)	SEQU	JENCI	E DES	CRIE	TION	l: SE	EQ II	NO:	:6:						
	Cys 1	Xaa	Xaa	His	Glu 5	Leu	Tyr	Val	Xaa	Phe 10	Xaa	Asp	Leu	Gly	Trp 15	Xaa
30	Asp	Trp	Xaa	Ile 20	Ala	Pro	Xaa	Gly	Tyr 25	Xaa	Ala	Tyr	Tyr	Cys 30	Glu	Gly
35	Glu	Cys	Xaa 35	Phe	Pro	Leu	Xaa	Ser 40	Xaa	Met	Asn	Ala	Thr 45	Asn	His	Ala
-	Ile	<b>Xaa</b> 50	Gln	Xaa	Leu	Val	His 55	Xaa	Xaa	Xaa	Pro	<b>X</b> aa 60	Xaa	Val	Pro	Lys
40	Xaa 65	Cys	Cys	Ala	Pro	Thr 70	Xaa	Leu	Xaa	Ala	Xaa 75	Ser	Val	Leu	Tyr	Xaa 80
	Asp	Xaa	Ser	Xaa	Asn 85	Val	Xaa	Leu	Xaa	Lys 90	Xaa	Arg	Asn	Met	Val 95	Val
45	Xaa	Ala	Cys	Gly 100	Cys	His										
	(2) INFOR	MATI	ON F	OR S	EQ I	D NO	:6:									
50	(i)			CHA												

PCT/US94/02335

10

- 54 -

- (B) TYPE: amino acid(C) STRANDEDNESS: single(D) TOPOLOGY: linear
- 5 (ii) MOLECULE TYPE: peptide
  - (ix) FEATURE:

    - (A) NAME/KEY: Cleavage-site
      (B) LOCATION: 1..4
      (D) OTHER INFORMATION: /note= "PROTEOLYTIC CLEAVAGE SITE"
  - (xi) SEQUENCE DESCRIPTION: SEQ ID NO:6:
- 15 Arg Xaa Xaa Arg

- 55 -

## What is claimed is:

1 1. An isolated binding partner having specific binding 2 affinity for an epitope on a soluble complex form 3 of an osteogenic protein.

4 5

6

7

8

9

10

said soluble complex form being characterized as a dimeric protein comprising a pair of polypeptide chain subunits associated to define a dimeric structure capable of inducing endochondral bone formation in a mammal when implanted in said mammal in association with a matrix, each said subunit having less than about 200 amino acids,

11 12 13

14

15

16

17 18

19

at least one of said subunits being noncovalently complexed with a peptide comprising a pro domain of a precursor form of an osteogenic protein subunit, or an allelic, species, or sequence variant thereof, to form a complex which is more soluble in aqueous solvents than the uncomplexed pair of subunits.

- 1 The isolated binding partner of claim 1 wherein
- 2 said binding partner is further characterized as
- 3 having substantially no binding affinity for the 4
  - mature, dimeric form of said osteogenic protein.
- 1 3. The isolated binding partner of claim 1 wherein
- said binding partner is further characterized as
- 3 having substantially no binding affinity for the
- isolated pro domain of one of said subunits. 4

- 56 -

An isolated binding partner having specific binding 1 4. affinity for an epitope on a soluble complex form 2 3 of an osteogenic protein. 4 said complex being characterized as a dimeric 6 protein comprising a pair of polypeptide chain 7 subunits associated to define a dimeric structure capable of inducing endochondral bone formation in 8 9 a mammal when implanted in said mammal in association with a matrix, each said subunit having 10 less than about 200 amino acids, 11 12 13 at least one of said subunits being 14 noncovalently complexed with a peptide comprising a 15 pro domain of a precursor form of an osteogenic 16 protein subunit, or an allelic, species, or sequence variant thereof, to form a complex which 17 is more soluble in aqueous solvents than the 18 19 uncomplexed pair of subunits, 20 21 said binding protein being characterized as 22 having affinity for said complex form and a said

pro domain of a precursor form, but not said

uncomplexed, dimeric protein.

23

24

- 57 -

An isolated binding partner having specific binding
 affinity for an epitope on a soluble complex form
 of an osteogenic protein,

said complex being characterized as a dimeric protein comprising a pair of protein subunits associated to define a dimeric structure capable of inducing endochondral bone formation in a mammal when implanted in said mammal in association with a matrix, each said subunit having less than about 200 amino acids.

at least one of said subunits being noncovalently complexed with a peptide comprising a pro domain of a precursor form of an osteogenic protein subunit, or an allelic, species, or sequence variant thereof, to form a complex which is more soluble in aqueous solvents than the uncomplexed pair of subunits.

said binding partner being characterized as having binding affinity for said complex form and said uncomplexed, dimeric form, but not a said pro domain of a precursor form of one of said subunits.

- An isolated binding partner having specific binding
   affinity for an epitope on a pro domain of an
- 3 osteogenic protein precursor form, said binding
- 4 partner having substantially no binding affinity
- for the mature, dimeric form or the soluble complex
- 6 form of said osteogenic protein.

- An isolated binding partner having specific binding
- affinity for an epitope on the mature dimeric form
- of an osteogenic protein, said binding partner 3
- 4 having substantially no binding affinity for the
- 5 soluble complex form or the pro domain of a
- precursor form of said osteogenic protein. 6
- The binding partner of claim 1, 4, 5, 6 or 7 1
- 2 wherein said binding partner is an antibody.
- 1 9. The antibody of claim 8 wherein said antibody is a
- 2 monoclonal antibody.
- 1 10. The antibody of claim 8 wherein said antibody is a
- 2 polyclonal antibody.
- 11. The binding partner of claim 1, 4, 5, 6 or 7 1
- wherein one said subunit of said dimeric osteogenic 2
- 3 protein is an OP-1 polypeptide chain including
- allelic, species and sequence variants thereof. 4
- 12. The binding partner of claim 11 wherein the other 1
- 2 subunit of said dimeric osteogenic protein is
- 3 selected from the group consisting of OP1, BMP2,
- BMP3 or BMP4, including allelic, species and 4
- 5 sequence variants thereof.
- 1 13. The binding partner of claim 11 wherein the other
- 2 subunit of said dimeric osteogenic protein is
- 3 selected from the group consisting of OP2, BMP5,
- 4 BMP6 or BMP9, including allelic, species and
- sequence variants thereof. 5

- 59 -

- 1 14. The binding partner of claim 11 wherein the other
- 2 subunit of said dimeric osteogenic protein is
- 3 selected from the group consisting of DPP, 60A,
- 4 Vgl, Vgr-1, including allelic, species and sequence
- 5 variants thereof.
- 1 15. The isolated binding partner of claim 1, 4, 5, 6
- 2 or 7 wherein said peptide comprises at least the
- 3 first 18 amino acids of an amino acid sequence
- 4 defining said pro domain.
- 1 16. The isolated binding partner of claim 15 wherein
- 2 said peptide comprises the full length form of said
- 3 pro domain.
- 1 17. The binding partner of claim 1, 4, 5, 6 or 7
- 2 wherein the pro domain peptide non-covalently
- 3 associated with said dimeric protein species
- 4 comprises an amino acid sequence selected from the
- 5 group consisting of OP1, OP2, BMP2, BMP3, BMP4,
- 6 BMP5, BMP6, DPP, Vg1, Vgr-1, 60A, including
- 7 allelic, species and sequence variants thereof.
- 1 18. The isolated binding partner of claim 1, 4, 5, 6 or
- 2 7 wherein said osteogenic protein is a synthetic
- 3 homolog of an osteogenic protein.
- 1 19. The isolated binding partner of claim 18 wherein
- 2 said synthetic homolog is selected from the group
- 3 consisting of COP 1, 3, 5, 7 and 16.
- 1 20. A method for specifically identifying a preselected
- 2 form of an osteogenic protein in solution, said
- 3 method comprising the steps of:

- 60 -

- 5 (a) exposing a solution likely to contain said
- 6 osteogenic protein to a binding partner having
- 7 specific binding affinity for a said preselected
- 8 form of a said osteogenically active protein under
- 9 conditions to promote specific binding between said
- 10 preselected protein form and said binding partner
- 11 to form a complex, and

12

- 13 (b)detecting the complex formed.
- 1 21. The method of claim 20 wherein said binding partner
- 2 is an antibody.
- 1 22. The method of claim 21 wherein said antibody is a
- 2 monoclonal antibody.
- 1 23. The method of claim 21 wherein said antibody is a
- 2 polyclonal antibody.
- 1 24. The method of claim 20 wherein said binding partner
- 2 has means for detection.
- 1 25. The method of claim 24 wherein said means for
- 2 detection comprises an enzyme or radioactive atom.
- 1 26. The method of claim 20 wherein said complex is
- 2 detected by means of a second binding partner
- 3 having specificity for said preselected osteogenic
- 4 protein.
- 1 27. The method of claim 20 wherein said preselected
- 2 form is present in admixture with other forms of
- 3 said protein.

- 61 -

- 28. The method of claim 20 comprising the additional
   step of quantitating the amount of complex formed.
- 1 29. The method of claim 20 wherein said solution in
- 2 step (a) is exposed to a first immobilized binding
- 3 partner having specificity for said preselected
- 4 form under conditions sufficient to promote
- 5 specific binding interaction between said first
- 6 binding partner and said preselected osteogenic
- 7 protein form to form a complex, and

8

- 9 said step of detecting (b) comprises the step of
- 10 exposing said complex to a second binding partner
- 11 having binding specificity for one of the proteins
- 12 in said complex.
  - 1 30. The method of claim 29 wherein said second binding
- 2 partner has specific binding affinity for said
- 3 preselected form of said osteogenic protein.
- 1 31. The method of claim 20 wherein said osteogenically
- 2 active protein is a dimeric protein comprising a
- 3 pair of disulfide bonded subunits, at least one of
- 4 said subunits comprising an OP-1 specific amino
- 5 acid sequence, including allelic, species and
- 6 sequence variants thereof.
- 1 32. The method of claim 31 wherein the other said
- 2 subunit comprises an amino acid sequence selected
- 3 from the group consisting of OP1, BMP2, BMP3 or
- 4 BMP4, including allelic, species and sequence
- 5 variants thereof.

- 62 -

1 2 3 4 5	33.	The method of claim 31 wherein the other said subunit comprises an amino acid sequence selected from the group consisting of OP2, BMP5, BMP6 or BMP9, including allelic, species and sequence variants thereof.
1	34.	The method of claim 31 wherein the other said
2		subunit comprises an amino acid sequence selected
3		from the group consisting of DPP, 60A, Vgl, Vgr-1,
4		including allelic, species and sequence variants
5		thereof.
1 2	35.	A kit for detecting a preselected form of an osteogenically active protein in solution, the kit
3		comprising:
4		
5		(a) means for capturing a fluid sample comprising
6		osteogenic protein,
7		
8		(b) a binding partner having specific binding
9		affinity for a said preselected form of said
10		osteogenic protein, and
11		
12		(c) means for detecting said binding partner bound
13		to said preselected osteogenic protein form.
1	36.	The kit of claim 35 wherein said binding partner is
2		an antibody.
1	27	The kit of claim 35 wherein said fluid sample is a
2	3/•	culture medium.
4		CATCAIC WEATON.

- 63 -

- 1 38. The kit of claim 35 wherein said means for
- 2 detecting (c) comprises a second binding partner
- 3 having specific binding affinity for said
- 4 preselected form of said osteogenic protein.
- 1 39. The kit according to claim 35 wherein said
- 2 osteogenically active protein is a dimeric protein
- 3 comprising a pair of disulfide bonded subunits, at
- 4 least one of said subunits comprising an OP-1
- 5 specific amino acid sequence, including allelic,
- 6 species and sequence variants thereof.
- 1 40. The kit according to claim 39 wherein said other
- 2 subunit is selected from the group consisting of
- 3 OP1, BMP2, BMP3 or BMP4, including allelic, species
- 4 and sequence variants thereof.
- 1 41. The kit according to claim 39 wherein said other
- 2 subunit is selected from the group consisting of
- 3 DPP, 60A, Vgl, Vgr-1, including allelic, species
- 4 and sequence variants thereof.
- 1 42. The kit according to claim 39 wherein said other
- 2 subunit is selected from the group consisting of
- 3 OP2, BMP5, BMP6 or BMP9, including allelic, species
- 4 and sequence variants thereof.





Fig. 1



BMP-3: RKKRSTGVLLPLQKSKNKKKQRKGPHRKSQTLQFDEQTLKKARRKQWIEPRNC	BMP-3: F
RCKRPRKRSYSKLPFTASNIC	Vg-1:
RISRSLPQGSGNWAQLRPLLVTFGHDGRGHALTRRRRRRSPKHHSQRARKKNKNC	BMP-4:
RHVRISRSLHQDEHSWSQIRPLLVTFGHDGRGHPLHRREKRQARHKQRKRLKSSC	BMP-2:
RSIRDVSGGEGGGKGGRNKRHARRPTRKNHDDTC	DPP:
RSKRSASHPRKRKKSVSPNNVPLLEPMESTRSC	60A:
RSVRAANKRKNQNRNKSSSHQDSSRMSSVGDYNTSEQKQAC	BMP-5:
RTTRSASSRRRQQSRNRSTQSQDVSRGSGSSDYNGSELKTAC	Vgr-1:
<u>RSIR</u> STGSKQRSQNRSKTPKNQEALRMANVAENSSSDQRQAC	OP-1:
RAPRSQOPFVVTFFRASPSPIRTPRAVRPLRRRQPKKSNELPQANRLPGIFDDVHGSHGRQVC	OP-2:

FIG.

Inte onal Application No PCT/US 94/02335

Relevant to claim No.

A. CLASSIFICATION OF SUBJECT MATTER IPC 5 C07K15/00 C12P21/08

C. DOCUMENTS CONSIDERED TO BE RELEVANT

G01N33/577

G01N33/68

According to International Patent Classification (IPC) or to both national classification and IPC

#### **B. FIELDS SEARCHED**

Minimum documentation searched (classification system followed by classification symbols)

IPC 5 CO7K C12P GO1N

Documentation searched other than minimum documentation to the extent that such documents are included in the fields searched

Electronic data base consulted during the international search (name of data base and, where practical, search terms used)

Citation of document, with indication, where appropriate, of the relevant passages

expression of a new member of the TGF-beta

see page 496, column 2 - page 497, column

superfamily in D. melanogaster'

1		
х	WO,A,93 00049 (ХОМА CORPORATION) 7 January 1993	1,3,5, 7-10,12, 17,18, 20-30,32
1	see page 57; example 8	
		ŀ
X	DEVELOPMENTAL BIOLOGY vol. 151 , 1992 pages 491 - 505	1-18, 20-34
	J. DOCTOR ET AL. 'Sequence, biochemical characterisation and developmental	

-/-

X	Further docume	nts are listed	in the contin	uation of box C.
* Spe	cial categories of	cited documen	nts :	

2, line 11 see page 501

- "A" document defining the general state of the art which is not considered to be of particular relevance
- "E" earlier document but published on or after the international filing date
- "L" document which may throw doubts on priority claim(s) or which is cited to establish the publication date of another citation or other special reason (as specified)
- "O" document referring to an oral disclosure, use, exhibition or other means
- "P" document published prior to the international filing date but later than the priority date claimed
- "T" later document published after the international filling date or priority date and not in conflict with the application but cited to understand the principle or theory underlying the invention

Patent family members are listed in annex.

- "X" document of particular relevance; the claimed invention cannot be considered novel or cannot be considered to involve an inventive step when the document is taken alone
- "Y" document of particular relevance; the claimed invention cannot be considered to involve an inventive step when the document is combined with one or more other such documents, such combination being obvious to a person skilled in the art.
- "&" document member of the same patent family

2 2. 07. 94

Date of the actual completion of the international search

Date of mailing of the international search report

12 July 1994

Authorized officer

Name and mailing address of the ISA

European Patent Office, P.B. 5818 Patentiaan 2
NL - 2220 HV Rijswijk
Tel. (- 31-70) 340-2040, Tx. 31 651 epo nl,
Fax: (+ 31-70) 340-3016

Skelly, J

Form PCT/ISA/218 (second sheet) (July 1992)

page 1 of 2

` 1

Inte onal Application No
PCT/US 94/02335

	PCT/US		
(Continue	anon) DOCUMENTS CONSIDERED TO BE RELEVANT  Citation of document, with indication, where appropriate, of the relevant passages	Relevant to claim No.	
acgory	Clearent of document, with minesuon, where appropriate, or the relevant passages		
	J. CELL. BIOL. vol. 120, no. 2 , 1993 pages 493 - 502 N. WALL ET AL. 'Biosynthesis and in vivo localisation of decapentaplegic Vg-related protein, DVR-6' see page 494, column 2 - page 495, column 2	6,8-10, 14	
	J. NEUROSCI. RES. vol. 33 , 1992 pages 142 - 147 U. SAUERMANN ET AL. 'Cloning of a novel TGF-beta related cytokine, the vgr, from rat brain' see page 142, column 2, line 30 - page 143, column 1, line 18	1,3,5, 7-10,14, 17,20	
(	PROC. NATL. ACAD. SCI. USA vol. 87 , 1990 pages 2220 - 2224 E. WANG ET AL. 'Recombinant huamn bone morphogenetic protein induces bone formation' see page 2220, column 2, line 15 - page 2222, column 1, line 15	1-18	
(	WO,A,89 09788 (CREATIVE BIOMOLECULES, INC) 19 October 1989 see page 31, line 10-20	7,18,19	
Ì	WO,A,92 15323 (CREATIVE BIOMOLECULES, INC.) 17 September 1992	1	
, P	WO,A,93 15197 (M.U.R.S.T.) 5 August 1993 see example 4	1-18	
(,P	WO,A,94 03600 (CREATIVE BIOMOLECULES, INC.) 17 February 1994 see the whole document	1-42	
	TRENDS IN GENETICS vol. 8, no. 3 , 1992 pages 97 - 102 V. ROSEN ET AL. 'The BMP proteins in bone formation and repair'	·	

1

aformation on patent family members

Inter vnal Application No PCT/US 94/02335

			. 01,700	
Patent document cited in search report	Publication date	Patent mem		Publication date
WO-A-9300049	07-01-93	US-A-	5284756	08-02-94
50000.5	0, 02 20	AŬ-A-	9162391	25-01-93
		CA-A-	2093790	21-12-92
		EP-A-	0546125	16-06-93
WO-A-8909788	19-10-89	US-A-	4968590	06-11-90
		US-A-	5002770	26-03-91
		US-A-	5011691	30-04-91
		AU-B-	628050	10-09-92
		AU-A-	3444989	03-11-89
		AU-B-	618357	19-12-91
		AU-A-	3530589	03-11-89
		EP-A-	0372031	13-06-90
		EP-A-	0362367	11-04-90
		JP-T-	3500655	14-02-91
		JP-T-	3502579	13-06-91
		WO-A-	8909787	19-10-89
		US-A-	5108753	28-04-92
		US-A-	5182365	26-01-93
		US-A-	5250302	05-10-93
		US-A-	5266683	30-11-93
		US-A-	5258494	02-11-93
		AU-B-	627850	03-09-92
		AU-A-	5174790	26-09-90
		EP-A-	0411105	06-02-91
		JP-T-	3504736	17-10-91
		WO-A-	9010018	07-09-90
		US-A-	4975526	04-12-90
		US-A-	5171574	15-12-92
		US-A-	5162114	10-11-92
WO-A-9215323			175 1000	06-10-92
	17-09-92	AUA-	1754392	
	17-09-92	CA-A-	2104678	12-09-92
	17-09-92			
 WO-A-9315197	17-09-92  05-08-93	CA-A-	2104678	12-09-92
WO-A-9315197 WO-A-9403600		CA-A- EP-A-	2104678	12-09-92
	05-08-93	CA-A- EP-A- NONE	2104678 0575555	12-09-92 29-12-93 

information on patent family members

Inter mal Application No
PCT/US 94/02335

Patent document ated in search report	Publication date				
WO-A-9403600		WO-A-	9403075	17-02-94	
		WO-A-	9403200	17-02-94	
		WO-A-	9406447	31-03-94	
		WO-A-	9406399	31-03-94	
		WO-A-	9406449	31-03-94	
		WO-A-	9406420	31-03-94	
		WO-A-	9410203	11-05-94	
•		AU-A-	3176293	27-04-93	
		CA-A-	2116559	01-04-93	
		EP-A-	0601135	15-06-94	
		WO-A-	9305751	01-04-93	

Form PCT/ISA/218 (patent family annex) (July 1992)